

# 24-Hour Emergency Number (ChemTel) 888.255.3924 in U.S.

## cleaning@alconox.com

## alconox.com

- distributors
- technical information
- free samples
- new developments
- white papers
- prices



# **LIQUINOX®** Critical-Cleaning Liquid Detergent

- Concentrated to save you money
- Replaces corrosive acids and hazardous solvents
- Phosphate-free, biodegradable and readily disposable
- Free rinsing to provide reliable results and no interfering residues
- For FDA cleaning validations and food contact surface cleaning

**Used to clean:** Healthcare instruments, laboratory ware, vacuum equipment, tissue culture ware, personal protective equipment, sampling apparatus, catheters, tubing, wine glasses, clean rooms, medical devices, optical parts, electronic components, pharmaceutical and cosmetics manufacturing equipment, metal castings, forgings and stampings, industrial parts, pipes, tanks and reactors. Used for phosphate sensitive analysis ware.

**Used to remove:** Soils, grit, grime, slime, grease, oils, blood, tissue, particulates, deposits, chemicals and solvents.

**Surfaces cleaned:** Corrosion inhibited formulation recommended for glass, metal, stainless steel, porcelain, ceramic, plastic, cement and fiber glass. Can be used on soft metals such as copper, aluminum, zinc and magnesium if rinsed promptly. Used for art restoration. Corrosion testing may be advisable.

**Cleaning method:** Soak, brush, sponge, cloth, ultrasonic, immersion/ recirculation cleaning in place. Will foam—not for spray or machine use.

**Directions:** Make a fresh 1% solution (2 1/2 tbsp per gal, 1 1/4 oz per gal, or 10 ml per L) in cold, warm or hot water. If available, use warm water. Use cold water for blood stains. For difficult soils, raise water temperature and use more detergent. Clean by soak, circulate,

wipe or ultrasonic method. Not for spray machines: will foam. RINSE THOROUGHLY—preferably with running water. For critical cleaning, do final or all rinsing in distilled, deionized or purified water. For food contact surfaces, rinse with potable water. Used on a wide range of glass, ceramic, plastic and metal surfaces. Corrosion testing may be advisable. Material compatibility guidelines available on request.

<b>Convenient Sizes</b>	Alcon	ox Cat. #			
Case 12 x 1 qt (12 x 0.95L)		1232		-	
Case of 4 x 1 gal (4 x 3	3.8L)	1201			
15 gal (57L) drum		1215		-	
55 gal (208L) drum		1255	-		
1 gal of concentrate makes 100 gal (1 L makes 100L) solution					

LIQUINOX passes inhibitory residue test for water analysis and is FDA certified. LIQUINOX is available from leading laboratory, hospital, clinical and industrial suppliers. To find a distributor for Alconox Inc. detergents, use the Locate Dealers search tool in the Support section of alconox.com. To request FREE samples, use the Get Samples form at alconox.com. Or call Alconox Inc. at +914.948.4040, today.

LQTB.2



## 24-Hour Emergency Number (ChemTel) 888.255.3924 in U.S.

## cleaning@alconox.com

#### alconox.com

- distributors
- technical information
- free samples
- new developments
- white papers
- prices

LIQUINOX Typical Data		<b>Typical Value</b>	
pH (as is)		8.5	
Surface Tension 1% Sol'n (dyn/cm)		29	
Conductivity 1% Sol'n	(mS/cm at 21 degrees C)	0.7	
Form		Liquid	
Color		Pale yellow	
Solubility in Water	Completely soluble in	all proportions	
Hard Water Effectivene	ess l	Highly effective	
Biodegradability	Yes, meet	s EC648/2004	
Foam Tendency		High foaming	
Shelf Life 24 r	months from end of month o	of manufacture	
Fragrance and Dye Cor	ntent	0%	
Phosphate Content (as phosphorus)		0%	
Flash Point (degrees F)		None	
Organic Carbon (% calculated w/w)		19%	
Density (lb/gal)		8.9	

## **Chemical Description**

Liquinox consists primarily of a homogeneous blend of sodium linear alkylaryl sulfonate, sodium xylene sulfonate, and ethoxylated alcohol. Liquinox contains a mixture of nonionic and anionic surfactants.

## **Cleaning Validation and Residue Detection**

To quantify trace levels of Liquinox for cleaning validation purposes, use conductivity, UV, TOC, or HPLC. For details see the Cleaning Validation page in Support at alconox.com. For lab accreditation tests to show the detergent is rinsed away, see the Accreditation Residue Detection in the Support section of alconox.com.

## LIQUINOX Technical Bulletin ©2022 Alconox Inc.

### Health and Safety Information

<b>OSHA Hazardous Ingre</b>	dients None
RCRA Hazard Class	Nonhazardous
<b>EPA Priority Pollutants</b>	None
DOT Hazard Class	Nonhazardous
Flammability	Nonflammable
Latex Content	None in product or packaging.
Irritation	Cat. 2 skin, Cat. 2A eye irritant
Oral Toxicity	$LD_{50} > 5000 \text{ mg/kg oral rat.}$
Inhalation Toxicity	Nonirritating
VOC Content	0%
GMO Content	None
BSE/TSE Source	None
2007/68/EC Annex IIIa	Allergens None
REACH Art. 57 (SVHC)	None
RoHS	No RoHS hazardous ingredients.
Carcinogenicity	NTP = NO IARC = NO OSHA = NO
Indgredients on TSCA I	nventory Yes
EC 907/2006 perfumes, all	Contains no enzymes, optical brighteners, ergenic fragrances, or preservative agents.

## **Precautions**

Avoid contact with skin, eyes and clothing. Do not take internally. In case of contact immediately flush skin with plenty of water. For eyes, flush with plenty of water for at least 15 minutes and get medical attention immediately. A Safety Data Sheet is available at alconox.com.

Contact Alconox Inc. for purchase specifications. Information presented is typical. Not to be taken as specifications. Typical data is not a specification.

While the information in this report should not be considered to be a product warranty, we urge you to investigate, test and verify the suitability of Alconox inc. detergents for your specific application. We, of course, can not give permission to use, or recommend the use of, our detergents where they infringe patents. No representation or warranty is made as to the safety of products or materials mentioned under the Federal Food Additives Amendment of 1958.

LQTB.2