

Distillation Units

Water Stills

Rivergum Series

Labglass

The Rivergum series of automatic all-glass Water Stills incorporates over a quarter of a century of Labglass experience in production of quality Water Stills for the world.

- Attractive Rivergum colour, powder coated finish
- 24 hour unattended operation of all models
- Can be supplied for use with pre-treated water sources



The Cascade Still produces 3 litre per hour single distilled water.

The Delta Still, unlike other stills available, allows re-configuration by the user at any time to produce 3 or 6 litres per hour single distilled water, or 3 litres per hour double distilled water.

Choosing a Water Still:

- Determine if double distilled water is required
If Yes - choose the Delta
- Determine production quantity of single distilled water required
If 6 litres per hour - choose the Delta
If 3 litres per hour - choose the Cascade



Rowe Code	Model	Output L/hr	Dual Feed Option	Double Distilled Possible
IS2160	Cascade	3	No	No
IS2165	Cascade	3	Yes	No
IS2170	Delta	3 or 6	No	Yes
IS2175	Delta	3 or 6	Yes	Yes