

Flasks

Boiling - Flat Bottom - narrow neck with rim.

Schott



Rowe Code	Capacity mL	Std. Pack Size
GF1345	100	10
GF1365	500	10
GF1435	2000	10
GF1440	3000	1

Boiling - Round Bottom - narrow neck with rim.

Schott




Rowe Code	Capacity mL	Std. Pack Size
GF14192	1000	10
GF14194	2000	10
GF1390	5000	1

Erlenmeyer - Narrow Neck

TGI

Borosilicate glass, graduated with rim



Rowe Code	Capacity mL	Neck ID mm	Base OD mm	Height mm	Std. Pack Size
GF1720	100	22	64	105	10
GF1825	200	34	79	131	10
GF1925	250	34	85	140	10
GF1985	300	34	87	156	10
GF2055	500	34	105	175	10
GF2135	1000	42	131	220	10
GF2220	2000	50	166	280	6
GF2425	5000	50	270	365	2

Erlenmeyer - Polypropylene

Kartell

For general laboratory use. Autoclavable.



Rowe Code	Capacity mL	Graduations mL	Neck NS	Std. Pack Size
PF1430	125	25	14/23	ea
PE1020	250	25	19/26	ea
PF1460	500	50	24/29	ea
PF1480	2000	100	34/35	ea

Contact us ...



ADELAIDE: 08 8186 0523 rowesa@rowe.com.au
 BRISBANE: 07 3376 9411 roweql@rowe.com.au
 MELBOURNE: 03 8795 7771 rowevic@rowe.com.au
 PERTH: 08 9302 1911 rowewa@rowe.com.au

www.rowe.com.au


Flasks

Erlenmeyer - Narrow Neck

Kimax

Borosilicate glass, graduated with reinforced beaded rim.

Rowe Code	Capacity mL	Base OD mm	Height mm	Std. Pack Size
GF1550	25	39	65	12
GF1650	50	50	78	12
GF1750	125	66	108	12
GF1900	250	81	130	12
GF1970	300	88	140	12
GF2100	500	102	174	6
GF2200	1000	128	213	6
GF2300	2000	161	275	4
GF2400	4000	208	355	1
GF2420	6000	235	393	1




Erlenmeyer - Wide Neck

TGI

Borosilicate glass, graduated with rim

Rowe Code	Capacity mL	Neck ID mm	Base OD mm	Height mm	Std. Pack Size
GF1600	50	34	51	85	10
GF2675	100	34	64	105	10
GF2815	250	50	85	135	10
GF2865	300	50	87	156	10
GF2905	500	50	105	170	10
GF2170	1000	50	131	220	10
GF2230	2000	72	153	280	6



Erlenmeyer - Wide Neck

Kimax

Borosilicate Glass, graduated with tooled top finish

Rowe Code	Capacity mL	Base OD mm	Height mm	Std. Pack Size
GF2651	125	66	108	12
GF2800	250	77	130	12
GF2950	500	97	172	6
GF2955	1000	122	216	6
GF2990	2000	161	265	4



Contact us ...



ADELAIDE: 08 8186 0523 rowesa@rowe.com.au
 BRISBANE: 07 3376 9411 roweqld@rowe.com.au
 MELBOURNE: 03 8795 7771 rowevic@rowe.com.au
 PERTH: 08 9302 1911 rowewa@rowe.com.au

www.rowe.com.au

Flasks




Filter (Buchner)

With tube connection. Borosilicate glass

TGI

Rowe Code	Capacity mL	Std. Pack Size
GF3265	500	ea
GF3345	1000	ea



Filter

Borosilicate Glass, Heavy wall, graduated with tube connection

Schott

Rowe Code	Capacity mL	Base OD mm	Height mm	Std. Pack Size
GF3105	100	65	100	10
GF3152	250	85	155	10
GF3260	500	105	185	10
GF3351	1000	135	230	10
GF3502	2000	166	265	1




Filter

Polypropylene. Moulded in one piece to withstand vacuum to 30in. Hg. Autoclavable. Tubulation at 45° angle prevents tipping.

Nalgene

Rowe Code	Capacity mL	Fits Tubing I.D.mm	Std. Pack Size
PF2203	500	6.4-9.5	ea
PF2205	1000	6.4-9.5	ea
PF2215	1700	9.5	ea



Kjeldahl

Rowe Code	Capacity mL	Joint Size	Std. Pack Size
GQ3165	250	24/29	ea
GQ3500	500	24/29	ea

Contact us ...



ADELAIDE: 08 8186 0523 rowesa@rowe.com.au
 BRISBANE: 07 3376 9411 roweqld@rowe.com.au
 MELBOURNE: 03 8795 7771 rowevic@rowe.com.au
 PERTH: 08 9302 1911 rowewa@rowe.com.au

www.rowe.com.au

Flask - Volumetric

Volumetric

Witeg

Borosilicate glass, Class A accuracy. Blue graduation. With polypropylene stopper.

Rowe Code	Capacity mL	Joint Size	Std. Pack Size
GF3660	5	10/19	ea
GF9020	10	10/19	ea
GF9002	20	10/19	ea
GF9003	25	10/19	ea
GF9004	50	12/21	ea
GF9006	100	12/21	ea
GF9005	200	14/23	ea
GF9007	250	14/23	ea
GF9008	500	19/26	ea
GF9009	1000	24/29	ea
GF9011	2000	29/32	ea
GF9012	5000	34/35	ea



Volumetric - Kohlrausch Pattern

Class A precision. Blue Graduation

Rowe Code	Capacity mL	Std. Pack Size
GF3621	100	ea




Volumetric - Polypropylene

Kartell

These flasks create no meniscus so final readings are easy. Calibrated gravimetrically at 20°C to Class B standards. Translucency of polypropylene enhanced when in contact with liquids. Thick walls eliminate distortion.

Thin necks make final reading precise. Autoclavable. Chemically non-absorbant. No wetting. Stoppers included.

Rowe Code	Capacity mL	OD mm	Neck Size	Std. Pack Size
PV1708	25	36	12/21	ea
PV1710	50	46	14/23	ea
PV1712	100	58	14/23	ea
PV1714	250	82	19/26	ea
PV1716	500	100	19/26	ea
PV1750	1000	120	19/26	ea



Contact us ...



ADELAIDE: 08 8186 0523 rowesa@rowe.com.au
 BRISBANE: 07 3376 9411 roweqld@rowe.com.au
 MELBOURNE: 03 8795 7771 rowevic@rowe.com.au
 PERTH: 08 9302 1911 rowewa@rowe.com.au

www.rowe.com.au