

pH Test Strips



Litmas paper



Universal paper



Duotest



pH - Fix

Litmus Paper

Rowe Code	Indicator Paper	Colour Change	pH Change	Std. Pack Size
FI1400	Litmus Paper Blue	Blue-Red	8.0 - 5.0	100 strips
FI1450	Litmus Paper Red	Red-Blue	5.0 - 8.0	100 strips

Macherey-Nagel

Universal Indicator Paper

Rowe Code	Indicator Paper	Std. Pack Size
FI1355	Universal Indicator pH 1-11 (1-2-3-4-5-6-7-8-9-10-11)	5m roll
FI1500	Universal Indicator pH 1-14 (1-2-3-5-6-7-8-9-10-12-14)	5m roll
FI1340	Special Indicator pH 6.4 - 8.0 (<6.4-6.4-6.6-6.8-7.0-7.2-7.4-7.6-7.8-8.0->8.0)	5m roll

Macherey-Nagel

Duotest indicator papers

Indicator papers with greater accuracy and increased ease of reading.

Rowe Code	Indicator Paper	Std. Pack Size
FI1700	pH 1-12 (1-2-3-4-5-6-7-8-9-10-11-12)	5m roll
FI1705	pH 3.5 - 6.8 (3.5-3.8-4.1-4.4-4.7-5.0-5.3-5.6-5.9-6.2-2.6-6.5-6.8)	5m roll
FI1710	pH 5.0 - 8.0 (5.0-5.3-5.6-5.9-6.2-6.5-6.8-7.1-7.4-7.7-8.0)	5m roll

Macherey-Nagel

Peñanon Indicator Papers Macherey-Nagel

for measurement of pH in coloured solutions

Rowe Code	Indicator Paper	Std. Pack Size
FI1065	pH 0.0 - 1.8 (0-0.3-0.6-0.8-1.0-1.2-1.5-1.8)	200 strips
FI1750	pH 10.5 - 13.0 (10.5-11.0-11.5-12.0-12.5-13.0)	200 strips

Note: other pH ranges available

pH - Fix Indicator Sticks,

Recommended for measurements of the pH of weakly buffered or strongly alkaline solutions. Non-bleeding

Macherey-Nagel

Rowe Code	Indicator Stick	Graduation of Scale	Std. Pack Size
FI1100	pH 1-14	0-1-2-3-4-5-6-7-8-9-10-11-12-13-14	100 sticks
FI1150	pH 0-6.0	0-0.5-1.0-1.5-2.0-2.5-3.0-3.5-4.0-4.5-5.0-5.5-6.0	100 sticks
FI1200	pH 2.0-9.0	2.0-2.5-3.0-3.5-4.0-4.5-5.0-5.5-6.0-6.5-7.0-7.5-8.0-8.5-9.0	100 sticks
FI1250	pH 4.5-10.0	4.5-5.0-5.5-6.0-6.5-7.0-7.5-8.0 8.5-9.0-9.5-10.0	100 sticks
FI1300	pH 7.0-14.0	7.0-7.5-8.0-8.5-9.0-9.5-10.0-10.5-11.0-11.5-12.0-12.5-13.0-13.5-14.0	100 sticks
FI1125	pH 0.3-2.3	0.3-0.7-1.0-1.3-1.6-1.9-2.3	100 sticks
FI1265	pH 5.1-7.2	5.1-5.4-5.7-6.0-6.3-6.6-6.9-7.2	100 sticks
FI1278	pH 6.0-7.7	6.0-6.4-6.7-7.0-7.3-7.7	100 sticks

Qualitative Test Papers

Test Papers allow the qualitative determination of ions and chemical compounds. They yield the information, whether a given component is present above the defined limit of detection.

Rowe Code	Indicator Paper	Std. Pack Size
FI1816	Chromium Test Paper (Chromium, Chromate)	200 strips
FI1825	Cyantesmo (Cyanides, Hydrocyanic Acid)	5m
FF6500	Fluoride Test Paper (Fluorides, Hydrogen Fluoride)	200 strips
FI1802	Dipyridyl Paper (Iron II)	200 strips
FI2100	Oil Test Paper (Oil in water & soil)	100 strips
FI1470	Peroxtesmo KO (Peroxidase in foodstuffs)	100 strips

Note: Other test papers available for a wide range of determinations.

pH Meter range
can be seen under
Meters

Contact us ...



ADELAIDE: 08 8186 0523 rowesa@rowe.com.au
 BRISBANE: 07 3376 9411 roweqld@rowe.com.au
 MELBOURNE: 03 8795 7771 rowevic@rowe.com.au
 PERTH: 08 9302 1911 rowewa@rowe.com.au

www.rowe.com.au