

## Riffle Splitters

We supply high quality Australian made as well as overseas sourced riffle splitters (Precision Carpco). They are available in mild steel and stainless steel and can be custom made to client specifications.

When ordering please specify:

- Mild steel or stainless steel
- Separation required eg  
50/50  
25/75  
12.5/87.5
- Riffle size, width and length
- Blade width
- Number of trays required
- Use in field or laboratory



Rowe Code	Description
RR1500	Riffle Splitter "Carpco" Stainless Steel 16x3mm Slots
RR1550	Riffle Splitter "Carpco" With Aluminium Trays 32x6mm Slots
RR1230	Riffle Splitter 6.35mm Slot Stainless Steel (Endecotts)

Contact us ...



ADELAIDE: 08 8186 0523 rowesa@rowe.com.au  
 BRISBANE: 07 3376 9411 roweqld@rowe.com.au  
 MELBOURNE: 03 8795 7771 rowevic@rowe.com.au  
 PERTH: 08 9302 1911 rowewa@rowe.com.au

[www.rowe.com.au](http://www.rowe.com.au)