

## Balances - Analytical



**Excellence Plus XP** analytical balances. Peak weighing performance and built-in QM-tools for full regulatory compliance



Item	Rowe Code	Description	Readability	Weighing Capacity
XP105DR	TBA	Balance, XP105 Delta Range®	0.01mg / 0.1mg	31g/120g
XP205	TBA	Balance, XP205	0.01mg	220g
XP205DR	IB3220	Balance, XP205 Delta Range®	0.01mg / 0.1mg	81g/220g
XP204	TBA	Balance, XP204	0.1mg	220g
XP504	TBA	Balance, XP504	0.1mg	520g

**Excellence XS** analytical balances. Innovative ergonomics for high weighing efficiency



XS105DU	IB3510	Balance, XS105DU Dual Range	0.01mg/0.1mg	41g/ 120g
XS205DU	IB3210	Balance, XS205DU Dual Range	0.01mg/0.1mg	81g/210g
XS64	TBA	Balance, XS64	0.1mg	61g
XS104	TBA	Balance, XS104	0.1mg	120g
XS204	IB3200	Balance, XS204	0.1mg	220g
XS204DR	TBA	Balance, XS204DR Delta Range	0.1mg/1mg	81g/ 220g

**Classic AB** analytical balances with FACT – for automatically accurate results



11139502	TBA	AB135-S/FACT Analytical Bal with FACT Dual range	0.01mg/ 0.1mg	31g/ 120g
11139506	TBA	AB265-S/FACT Analytical Bal with FACT Dual range	0.01mg/ 0.1mg	61g/ 220g
11139574	TBA	AB54-S/FACT Analytical Balance with FACT	0.1mg	51g
11139578	TBA	AB104-S/FACT Analytical Balance with FACT	0.1mg	110g
11139582	IB3520	AB204-S/FACT Analytical Balance with FACT	0.1mg	220g
11139586	TBA	AB304-S/FACT Analytical Balance with FACT	0.1mg	320g

**Classic Light AL** analytical balance. Low-budget balance for simple weighing tasks



AL54/00	IB3505	Balance	0.1 mg	51 g
AL54/01	TBA	Balance with RS232	0.1 mg	51 g
AL104/00	TBA	Balance	0.1 mg	110 g
AL104/01	TBA	Balance with RS232	0.1 mg	110 g
AL204/00	TBA	Balance	0.1 mg	210 g
AL204/01	TBA	Balance with RS232	0.1 mg	210 g

Optional RS232 interface supplied with external calibration weight.

Contact us ...



ADELAIDE: 08 8186 0523 rowesa@rowe.com.au  
 BRISBANE: 07 3376 9411 roweql@rowe.com.au  
 MELBOURNE: 03 8795 7771 rowevic@rowe.com.au  
 PERTH: 08 9302 1911 rowewa@rowe.com.au

[www.rowe.com.au](http://www.rowe.com.au)

## Balances - Analytical



**Analytical Plus** range of durable, reliable balances delivers performance and precision you would expect to find in premium priced balances.



Item	Rowe Code	Readability	Weighing Capacity
AP250D	IB7050	0.01 mg / 0.1 mg	52 g / 210 g
AP110	TBA	0.1 mg	110 g
AP210	TBA	0.1 mg	210 g
AP310	TBA	0.1 mg	310 g

**Voyager Pro's** advanced applications simplify complex laboratory measurements and help reduce errors in manual calculations.



EP64	TBA	0.1 mg	62 g
EP64C	TBA	0.1 mg	62 g
EP214D	TBA	0.1 mg / 1 mg	100 g / 210 g
EP214DC	TBA	0.1 mg / 1 mg	100 g / 210 g
EP114	TBA	0.1 mg	110 g
EP114C	TBA	0.1 mg	110 g
EP214	IB7414	0.1 mg	210 g
EP214C	IB7415	0.1 mg	210 g

**Adventurer Pro** Offers more standard features and options than any other balance in its class, the Adventurer Pro has the flexibility to meet your application needs.



AV114	TBA	0.1 mg	110 g
AV114C	TBA	0.1 mg	110 g
AV264	TBA	0.1 mg	260 g
AV264C	IB7146	0.1 mg	260 g

**Adventurer SL.** Designed for uncomplicated performance, the Adventurer SL has only those features necessary for routine weighing.



AS214	TBA	0.1 mg	210 g
-------	-----	--------	-------

**Adventurer. Ohaus®** Analytical balances offer the right range of capacities and functions to meet your most demanding requirements.



AR0640	TBA	0.1 mg	65 g
AR1140	TBA	0.1 mg	110 g
AR2140	IB7145	0.1 mg	210 g

Contact us ...



ADELAIDE: 08 8186 0523 rowesa@rowe.com.au  
 BRISBANE: 07 3376 9411 roweqld@rowe.com.au  
 MELBOURNE: 03 8795 7771 rowevic@rowe.com.au  
 PERTH: 08 9302 1911 rowewa@rowe.com.au

[www.rowe.com.au](http://www.rowe.com.au)

## Balances - Bench Scales



### CKW Series

The CKW Series all stainless steel construction scale is designed for check weighing applications requiring fast and accurate operation. Combining a NSF/USDA certified food-safe design, durable stainless steel construction, and industry leading software, the Ohaus CKW series of bench scales are dedicated check weighing scales built for applications requiring fast and accurate operation. Designed with the operator in mind, the Ohaus CKW scales are easy to operate, but allow full configuration to match different checkweighing procedures.



Item	Rowe Code	Readability	Weighing Capacity
CKW3R55	TBA	0.0005 kg 0.001 lb	3 kg 6 lb
CKW6R55	TBA	0.001 kg 0.002 lb	6 kg 15 lb
CKW15L55	TBA	0.002 kg 0.005 lb	15 kg 30 lb
CKW30L55	TBA	0.005 kg 0.01 lb	30 kg 60 lb
<b>CKW Indicator</b>			
CKW55	TBA	1:10,000 (NTEP pending)	1:20,000 max
<b>CKW Bases</b>			
CKW3R	TBA	0.0005 kg	3 kg
	TBA	0.001 lb	6 lb
CKW6R	TBA	0.001 kg	6 kg
	TBA	0.002 lb	15 lb
CKW15L	TBA	0.002 kg	15 kg
	TBA	0.005 lb	30 lb
CKW30L	TBA	0.005 kg	30 kg

Contact us ...



ADELAIDE:  
BRISBANE:  
MELBOURNE:  
PERTH:

08 8186 0523  
07 3376 9411  
03 8795 7771  
08 9302 1911

rowesa@rowe.com.au  
roweqld@rowe.com.au  
rowevic@rowe.com.au  
rowewa@rowe.com.au

[www.rowe.com.au](http://www.rowe.com.au)

## Balances - Bench Scales

Item	Rowe Code	Readability	Weighing Capacity
ES Series			
ES6R	TBA	0.002 kg	6 kg
ES30R	IB7590	0.01 kg	30 kg
ES50L	TBA	0.02 kg	50 kg
ES50R	TBA	0.02 kg	50 kg
ES100L	TBA	0.05 kg	100 kg
ES200L	TBA	0.1 kg	200 kg



Utility and economy come together in the ES Series of all-purpose bench scales. The low profile design of the ES Series scales makes them perfect for shipping, receiving, and general bench applications.



### Trooper

The Trooper's design is optimized for economical weighing applications. A concave stainless steel weighing pan and ABS housing provide protection for light industrial conditions, and the simple menu operation and easy-to-read backlit LCD display make the Trooper an affordable weighing solution. Portability and compact footprint allow use of the Trooper throughout your business. (60lbs) Readability: 0.0005kg to 0.005kg

Item	Rowe Code	Readability	Weighing Capacity
TR3RS	TBA	0.0005 kg	3 kg
TR6RS	IB3723	0.001 kg	6 kg
TR15RS	TBA	0.002 kg	15 kg
TR30RS	IB7585	0.005 kg	30 kg



Ohaus Corporation is a leading manufacturer of scales and balances for the laboratory, education, industrial and speciality markets worldwide.

With over 90 years of experience in the design, development and marketing of balances for a wide variety of applications, Ohaus has earned its reputation as a trusted supplier of precise, reliable and affordable products that meet virtually any weighing need. We remain firmly focused on continuing the long-standing record of success and on expanding our global presence with new products, new solutions and new relationships worldwide.

Contact us ...



ADELAIDE: 08 8186 0523 rowesa@rowe.com.au  
 BRISBANE: 07 3376 9411 roweqld@rowe.com.au  
 MELBOURNE: 03 8795 7771 rowevic@rowe.com.au  
 PERTH: 08 9302 1911 rowewa@rowe.com.au

[www.rowe.com.au](http://www.rowe.com.au)

## Balances - Bench Scales



### Champ II

Ohaus® Champ II Bench Scales are both durable and economical. These portable, 'plug-and-play' scales are ready to go to work for you right out of the box.



Item	Rowe Code	Readability	Weighing Capacity
CH15R11	TBA	0.002 kg	15 kg
CH30R11	IB7039	0.005 kg	30 kg
CH60L11	IB7070	0.01 kg	60 kg
CH60R11	TBA	0.01 kg	60 kg
CH100R11	TBA	0.02 kg	100 kg
CH150R11	TBA	0.02 kg	150 kg
CH300R11	TBA	0.05 kg	300 kg

### Champ General

For general weighing needs, there's a Champ Bench Scale that's just right for your application and budget.



CQ10R11	TBA	0.001 kg	10 kg
CQ25R11	TBA	0.002 kg	25 kg
CQ50L11	IB7038	0.005 kg	50 kg
CQ100L11	TBA	0.01 kg	100 kg
CQ250XL11	TBA	0.02 kg	250 kg

### Champ Multi-functional

From multiple weighing needs such as shipping to precise industrial requirements, there's a Champ Multi-Functional Bench Scale that's just right for your application and budget.



CQ10R31	TBA	0.001 kg	10 kg
CQ10R33	TBA	0.001 kg	10 kg
CQ25R31	IB7010	0.002 kg	25 kg
CQ25R33	TBA	0.002 kg	25 kg
CQ50L31	TBA	0.005 kg	50 kg
CQ50L33	TBA	0.005 kg	50 kg
CQ100L31	TBA	0.01 kg	100 kg
CQ100L33	TBA	0.01 kg	100 kg
CQ250XL31	TBA	0.02 kg	250 kg
CQ250XL33	TBA	0.02 kg	250 kg

### Champ Washdown

From general weighing needs such as shipping to precise industrial and water washdown requirements, there's a Champ Bench Scale that's just right for your application and budget.



CQ10R11W	IB7037	0.001 kg	10 kg
CQ25R11W	TBA	0.002 kg	25 kg
CQ50L11W	IB7036	0.005 kg	50 kg
CQ100L11W	TBA	0.01 kg	100 kg
CQ250XL11W	TBA	0.02 kg	250 kg

Contact us ...



ADELAIDE:  
BRISBANE:  
MELBOURNE:  
PERTH:

08 8186 0523  
07 3376 9411  
03 8795 7771  
08 9302 1911

rowesa@rowe.com.au  
roweqld@rowe.com.au  
rowevic@rowe.com.au  
rowewa@rowe.com.au

[www.rowe.com.au](http://www.rowe.com.au)

## Balances - Counting



Item	Rowe Code	Readability	Weighing Capacity
EC3	IB7060	0.1 g	3000 g
EC6	TBA	0.2 g	6000 g
EC15	TBA	0.5 g	15000 g
EC30	IB7065	1 g	30000 g



The **EC Series** is designed for fast, accurate operations ideally suited for annual inventories, parts rooms, rentals and other counting applications. Easy to use, the Ohaus® EC Series offers accuracy, durability and portability in one practical design.

TC3RS	TBA	0.0005 kg	3 kg
TC6RS	TBA	0.001 kg	6 kg
TC15RS	TBA	0.002 kg	15 kg
TC30RS	TBA	0.005 kg	30 kg



The **Trooper Count** is an economical counting scale with enhanced software designed for flexible counting operations. Featuring simple to advanced counting modes, 300 memory locations for up to 60 items, and an intuitive, easy-to-operate design, the Ohaus® Trooper Count is a practical counting product you can depend on for a full range of industrial use.

RC3RS	TBA	0.0002 kg	3 kg
RC6RS	TBA	0.0005 kg	6 kg
RC12LS	TBA	0.001 kg	12 kg
RC30LS	TBA	0.002 kg	30 kg
RC60LS	TBA	0.005 kg	60 kg



The **Ranger Count** series is a comprehensive line of high performance counting bench scales.

## Balances - Mechanical



### Triple Beam

Whether you're weighing solids, liquids, powders or even animals, an Ohaus® Triple Beam balance is well-equipped to handle the job.



### Harvard Trip

With a two-pan design, the Harvard Trip balances feature magnetic damping, beams with sliding weights and a variety of weighing platforms



**Dial-O-Gram® & Cent-O-Gram®** Favorites in classrooms and industries worldwide, these balances offer high precision and remarkable value.



### Field Test Scale

Ohaus® **Field Test** Scales have set an industry standard for precision in high-capacity weighing. Durable enough to resist corrosion, these precise balances also weigh with a feather touch - even under the most demanding conditions



### Ohaus® Heavy Duty

Solution Balance has set an industry standard for precision in high-capacity weighing.



Item	Rowe Code	Readability	Weighing Capacity
1610-00	IB1706	0.1 g / 0.1 g	610 g / 2610 g
1650-00	TBA	0.1 g / 0.1 g	610 g / 2610 g
710-00	TBA	0.1 g / 0.1 g	610 g / 2610 g
710-T0	TBA	0.1 g / 0.1 g	610 g / 2610 g
720-00	TBA	0.1 g / 0.1 g	610 g / 2610 g
720-S0	TBA	0.1 g / 0.1 g	610 g / 2610 g
730-00	TBA	0.1 g / 0.1 g	610 g / 2610 g
750-S0	TBA	0.1 g / 0.1 g	610 g / 2610 g
760-00	TBA	0.1 g / 0.1 g	610 g / 2610 g
1650-W0	IB1705	0.1 g	2610 g
750-SW	IB1710	0.1 g	2610 g

1450-SD	IB1730	0.1 g	2000 g
1550-SD	TBA	0.1 g	2000 g
1560-SD	TBA	0.1 g	2000 g

310-00	TBA	0.01 g	310 g
311-00	TBA	0.01 g	311 g

2400-11	TBA	5 g	16 kg
2400-12	TBA	5 g	16 kg

Ohaus® Field Test Scales have set an industry standard for precision in high-capacity weighing. Durable enough to resist corrosion, these precise balances also weigh with a feather touch - even under the most demanding conditions

1119-D0	TBA	1 g	20 kg
---------	-----	-----	-------

Durable enough to resist corrosion, this precise balance also weighs with a feather touch - even under the most demanding conditions.

Contact us ...



ADELAIDE:  
BRISBANE:  
MELBOURNE:  
PERTH:

08 8186 0523  
07 3376 9411  
03 8795 7771  
08 9302 1911

rowesa@rowe.com.au  
roweqld@rowe.com.au  
rowevic@rowe.com.au  
rowewa@rowe.com.au

[www.rowe.com.au](http://www.rowe.com.au)

## Balances - Micro



### Microbalances

The METTLER TOLEDO UMX is an innovative micro balance, which sets new world standards for weighing performance, user friendliness and quality standards. The UMX has a very high resolution: 50 million digits. A brilliant touch screen display, intelligent user guidance, and hands-free infrared draft shield control make this highly-precise balance faster, easier and more fun to use than ever before

Item	Rowe Code	Description	Readability	Weighing Capacity
UMX2	TBA	Microbalance, UMX2	0.1µg	2100mg
MX5	IB2092	Microbalance, MX5	1µg	5100mg
AX26DR	TBA	Balance, AX26 Delta Range®	2µg / 0.01mg	3g/21g



Contact us ...



ADELAIDE: 08 8186 0523 rowesa@rowe.com.au  
 BRISBANE: 07 3376 9411 roweql@rowe.com.au  
 MELBOURNE: 03 8795 7771 rowevic@rowe.com.au  
 PERTH: 08 9302 1911 rowewa@rowe.com.au

[www.rowe.com.au](http://www.rowe.com.au)

## Balances - Moisture Analysis



The halogen moisture analyzers from METTLER TOLEDO warrant reliable moisture determination within just minutes. The circular halogen lamp and gold-plated reflector ensure the exceptionally fast and repeatable heating of the sample. The easy operation, automatic method development and sophisticated QM Tools are among the unique features for regulated use.



Item	Rowe Code	Description	Readability	"Weighing Capacity"
HR83P	TBA	Halogen Moisture analyser, HR83, with built in printer	1mg	81g
HR83	IB2038	Halogen Moisture analyser, HG83	1mg	81g
HG63P	IB2040	Halogen Moisture analyser, HG63, with built in printer	1mg	61g
HG63	TBA	Halogen Moisture analyser, HG63	1mg	61g
HB43	IB2012	Halogen Moisture analyser, HB43	1mg	41g



Contact us ...



ADELAIDE: 08 8186 0523 rowesa@rowe.com.au  
 BRISBANE: 07 3376 9411 roweqld@rowe.com.au  
 MELBOURNE: 03 8795 7771 rowevic@rowe.com.au  
 PERTH: 08 9302 1911 rowewa@rowe.com.au

[www.rowe.com.au](http://www.rowe.com.au)

## Balances - Moisture Analysis

Designed to quickly, effectively and affordably measure moisture content, the Ohaus MB Series Moisture Analyzers incorporate state-of-the-art halogen heating with precision weighing technology to give a fast and accurate method for moisture content determination.



Perfect for applications in the pharmaceutical, chemical and research industries; versatile and rugged enough for continuous operation in food and beverage, quality control, environmental and many other applications.

Item	Rowe Code	Readability	Weighing Capacity
MB45	IB2019	0.001g	45g
MB35	IB2018	0.005g	35g



Contact us ...



ADELAIDE: 08 8186 0523 rowesa@rowe.com.au  
 BRISBANE: 07 3376 9411 roweqld@rowe.com.au  
 MELBOURNE: 03 8795 7771 rowevic@rowe.com.au  
 PERTH: 08 9302 1911 rowewa@rowe.com.au

[www.rowe.com.au](http://www.rowe.com.au)

## Balances - Moisture Analysis Accessories

### Accessories - Moisture Analyser

Sample Handling		
214526	TBA	Pan handler, sample, set of 3
13865	RT3555	Aluminium dishes, box of 80 pcs for one time use
11113863	TBA	Aluminium dishes professional thick walled box of 80 pcs
214462	TBA	Pans, HA-DR1, reusable, stainless steel, set of 3 (6mm high)
13954	TBA	Pans, reusable, stainless steel , set of 3 (15mm high)
214695	TBA	Cage for high volume lightweight material
214464	TBA	Fibreglass paper, box of 100 pcs
Printers		
HA-P43	TBA	Printer, HA-P43 (built-in) for HR/HG only
RS-P42	TBA	Printer, RS-P42
65975	TBA	Ink ribbon, black
72456	TBA	Paper, set of 5 rolls
11600388	TBA	Paper, self-adhesive, 1 roll
Quality Management		
214455	TBA	Temperature calibration/adjustment kit, HA-TC
214528	TBA	Temperature calibration/adjustment kit, certified, HA-TCC
214534	TBA	Recalibration, HA-TCC (inc. retraceable test cert. at 50C, 100C, 160C)
00158640	TBA	Single weight class F1, 20 g for HB
00158650	TBA	Single weight class F1, 50 g for HR/HG
Accessories		
21900879	TBA	Barcode scanner RS232 (must include power supply and cable)
214533	TBA	In-use cover, set of 2 for HR/HG
214515	TBA	Carry case for HR/HG
11113363	TBA	In-use cover, set of 2 for HB
11113855	TBA	Carry case for HB
11600361	TBA	Antitheft device (steel cord with lock) for HB
214469	TBA	Dryer replacement unit
13701	TBA	Thermometer kit, LP16 / LJ16
214240	TBA	Temperature calibration kit, LP16 / LJ16
214247	TBA	Thermal strips for 00214240
237010	TBA	BalanceLink Software 3.5" Disk (windows and NT)
21900855	TBA	BalanceLink Upgrade Version 2.2 to Version 3.0
Cables		
229050	TBA	Cable, LC-RS25,RS232/25 Pin, m/f, for IBM PC and compatibles 1.5 m
229065	TBA	Cable, LC-RS9, RS232/9 Pin, m/f, for IBM PC and compatibles, 1.5 m
11101051	TBA	Cable, RS9m – RS9f , 2 metre for HB
11101052	TBA	Cable, RS9m – RS25f , 2 metre for HB

Contact us ...



ADELAIDE: 08 8186 0523 rowesa@rowe.com.au  
 BRISBANE: 07 3376 9411 roweqld@rowe.com.au  
 MELBOURNE: 03 8795 7771 rowevic@rowe.com.au  
 PERTH: 08 9302 1911 rowewa@rowe.com.au

[www.rowe.com.au](http://www.rowe.com.au)

## Balances - Portable



Item	Rowe Code	Readability	Weighing Capacity
SP202	IB7710	0.01 g	200 g
SP401	IB7705	0.1 g	400 g
SP402	IB7720	0.01 g	400 g
SP601	TBA	0.1 g	600 g
SP2001	IB7760	0.1 g	2000 g
SP4001	IB7770	0.1 g	4000 g
SP6000	IB7765	1 g	6000 g

The **Scout Pro** combines the accuracy and reliability of earlier Scout balances with all new features for an exceptional balance at an affordable price.



AV53	IB7045	0.001 g	51 g
AV212	TBA	0.01 g	210 g
AV412	IB7182	0.01 g	410 g
AV812	IB7185	0.01 g	810 g
AV2101	TBA	0.1 g	2100 g
AV4101	IB7190	0.1 g	4100 g
AV8101	IB7580	0.1 g	8100 g

Ohaus **Adventurer Pro** balances are ideal for general laboratory, industrial and educational applications requiring multiple weighing units. Ohaus Adventurer Pro balances are the economical way to weigh with precision.



AS153	TBA	0.001 g	150 g
AS312	TBA	0.01 g	310 g
AS612	TBA	0.01 g	610 g
AS811	TBA	0.1 g	810 g
AS1502	TBA	0.01 g	1500 g
AS3101	TBA	0.1 g	3100 g
AS6101	TBA	0.1 g	6100 g
AS8100	TBA	1 g	8100 g

Designed for uncomplicated performance, the **Adventurer SL** has only those features necessary for routine weighing.



Contact us ...



ADELAIDE:	08 8186 0523	rowesa@rowe.com.au
BRISBANE:	07 3376 9411	roweqld@rowe.com.au
MELBOURNE:	03 8795 7771	rowevic@rowe.com.au
PERTH:	08 9302 1911	rowewa@rowe.com.au

[www.rowe.com.au](http://www.rowe.com.au)

## Balances - Portable



The Ohaus® Navigator balance is a full-featured, compact, affordable, portable balance. The Navigator balance offers you capacities from 32g to 8100g.

(Note: Balances beginning with item # N1xxxx have Internal calibration and those with N0xxxx do not.)

Item	Rowe Code	Readability	Weighing Capacity
N10330	TBA	0.002 g	32 g
N12120	TBA	0.01 g	210 g
N14120	TBA	0.01 g	410 g
N06112	TBA	0.1 g	610 g
N18110	TBA	0.1 g	810 g
N1B110	IB7225	0.1 g	2100 g
N1D110	TBA	0.1 g	4100 g
N0F1A2	TBA	1 g	6100 g
N1H110	TBA	0.5 g	8100 g

Ohaus® **Compact Scales** are the perfect economical portable balance. Ohaus® Compact Scales are ideal when cost matters and dependable results are a must, offering capacities of 200g to 5000g and readabilities from 0.1g. "P" denotes US Postal version.



CS200	IB1810	0.1 g	200 g
CS2000	IB1840	1 g	2000 g
CS2000P	TBA	1 g	2000 g
CS5000	TBA	2 g	5000 g
CS5000P	TBA	2 g	5000 g

Ohaus® **Hand-Held Scales** are convenient and compact allowing weighing in the palm of your hand. Use them to weigh jewellery, take them out into the field, for quality control, these scales are ideal to carry anywhere. They easily fit into your pocket or backpack for weighing on the go with capacities up to 320g and readabilities of 0.1g.



HH120D	TBA	0.1 g / 0.2 g	60g / 120 g
HH120	TBA	0.1 g	120 g
HH320	IB7180	0.1 g	320 g

**Pocket Scales** are the smallest and most portable electronic scales from Ohaus®. Rugged and innovative design lets you take the precision of an Ohaus® scale anywhere with capacities up to 250g and readabilities of 0.1g.



PS121	TBA	0.1 g	120 g
PS251	TBA	0.1 g	250 g

Contact us ...



ADELAIDE: 08 8186 0523  
BRISBANE: 07 3376 9411  
MELBOURNE: 03 8795 7771  
PERTH: 08 9302 1911

rowesa@rowe.com.au  
roweqld@rowe.com.au  
rowevic@rowe.com.au  
rowewa@rowe.com.au

[www.rowe.com.au](http://www.rowe.com.au)

## Balances - Precision


**Excellence precision balances XP S platforms**

Leading in terms of weighing performance and ease of use, our Excellence balances increase your productivity, support regulatory compliance and offer clever solutions for instrument control and data management.



Item	Rowe Code	Description	Readability	Weighing Capacity
XP204S	TBA	Balance, XP204S	0.1mg	210g
XP404S	TBA	Balance, XP404S	0.1mg	410g
XP404SDR	TBA	Balance, XP404S Delta Range®	0.1mg/1mg	80g/410g
XP203S	TBA	Balance, XP203S	0.001g	210g
XP603S	TBA	Balance, XP603S	0.001g	610g
XP603SDR	TBA	Balance, XP603S Delta Range®	0.001g / 0.01g	120g/ 610g
XP1203S	TBA	Balance, XP1203S	0.001g	1210g
XP2003SDR	IB3710	Balance, XP2003S Delta Range®	0.001g / 0.01g	500g/ 2100g
XP5003SDR	TBA	Balance, XP5003S Delta Range®	0.001g / 0.01g	1000g/ 5100g
XP1202S	TBA	Balance, XP1202S	0.01g	1210g
XP4002S	TBA	Balance, XP4002S	0.01g	4100g
XP6002S	TBA	Balance, XP6002S	0.01g	6100g
XP6002SDR	TBA	Balance, XP6002S Delta Range®	0.01g / 0.1g	1200g/ 6100g
XP8002S	TBA	Balance, XP8002S	0.01g	8100g
XP10002S	TBA	Balance, XP10002S	0.01g	10100g
XP10002SDR	TBA	Balance, XP10002S Delta Range®	0.01g / 0.1g	2000g/ 10100g
XP4001S	TBA	Balance, XP4001S	0.1g	4100g
XP6001S	TBA	Balance, XP6001S	0.1g	6100g
XP8001S	TBA	Balance, XP8001S	0.1g	8100g
XP10001S	TBA	Balance, XP10001S	0.1g	10100g

**Excellence precision balances XP M platforms**


XP6002MDR	TBA	Balance, XP6002M Delta Range®	0.01g / 0.1g	1200g/ 6100g
XP12002MDR	TBA	Balance, XP12002M Delta Range®	0.01g / 0.1g	2400g/ 12100g
XP8001M	TBA	Balance, XP8001M	0.1g	8100g
XP8001MDR	TBA	Balance, XP8001M Delta Range®	0.1g / 1.0g	1600g/ 8100g
XP12001M	TBA	Balance, XP12001M	0.1g	12100g
XP12000M	TBA	Balance, XP12000M	1.0g	12100g
XP16001M	TBA	Balance, XP16001M	0.1g	16,100g
XP20001M		Balance, XP20001M	0.1g	20,100g
XP20000M		Balance, XP20000M	1g	20,100g

Contact us ...



ADELAIDE: 08 8186 0523 rowesa@rowe.com.au  
 BRISBANE: 07 3376 9411 roweqld@rowe.com.au  
 MELBOURNE: 03 8795 7771 rowevic@rowe.com.au  
 PERTH: 08 9302 1911 rowewa@rowe.com.au

[www.rowe.com.au](http://www.rowe.com.au)

## Balances - Precision



## Excellence precision balances XP L platforms



Item	Rowe Code	Description	Readability	Weighing Capacity
XP8001L	TBA	Balance, XP8001	0.1g	8,100g
XP16001L	TBA	Balance, XP16001L	0.1g	16,100g
XP32001L	TBA	Balance, XP32001	0.1g	32,100g
XP32001LDR	TBA	Balance, XP32001L Delta Range®	0.1g / 1g	6400g/ 32100g
XP64001L	TBA	Balance, XP64001L	0.1g	64,100g
XP16000L	TBA	Balance, XP16000L	1g	16,100g
XP32000L	TBA	Balance, XP32000L	1g	32,100g
XP64000L	TBA	Balance, XP64000L	1g	64100g

## Excellence precision balances XS S platforms



XS203S	TBA	Balance, XS203S	0.001g	210 g
XS403S	TBA	Balance, XS403S	0.001g	410 g
XS603S	TBA	Balance, XS603S	0.001g	610 g
XS603SDR	TBA	Balance, XS603S DeltaRange®	0.001g / 0.01g	120 g / 610 g
XS1003S	TBA	Balance, XS1003S	0.001g	1010 g
XS802S	TBA	Balance, XS802S	0.01g	810 g
XS2002S	TBA	Balance, XS2002S	0.01g	2100 g
XS4002S	TBA	Balance, XS4002S	0.01g	4100 g
XS4002SDR	TBA	Balance, XS4002S DeltaRange®	0.01g / 0.1g	800g/4100 g
XS6002S	TBA	Balance, XS6002S	0.01g	6100 g
XS6002SDR	TBA	Balance, XS6002S DeltaRange®	0.01g / 0.1g	1200g/ 6100g
XS4001S	TBA	Balance, XS4001S	0.1g	4100g
XS6001S	TBA	Balance, XS6001S	0.1g	6100g
XS8001S	TBA	Balance, XS8001S	0.1g	8100g

## Excellence precision balances XS M platforms



XS6001M	TBA	Balance, XS6001M	0.1g	6100g
XS6001MDR	TBA	Balance, XS6001M DeltaRange®	0.1g / 1g	1200g/ 6100g
XS10001M	TBA	Balance, XS10001M	0.1g	10100g
XS10000M	TBA	Balance, XS10000M	1g	10100g
XS16001M	TBA	Balance, XS16001M	0.1g	16,100g
XS16000M	TBA	Balance, XS16000M	1g	16,100g

## Excellence precision balances XS L platforms



XS8001L	TBA	Balance, XS8001	0.1g	8,100g
XS16001L	TBA	Balance, XS16001L	0.1g	16,100g
XS32001L	TBA	Balance, XS32001	0.1g	32,100g
XS32001LDR	TBA	Balance, XS32001L Delta Range®	0.1g / 1g	6400g/ 32100g
XS16000L	TBA	Balance, XS16000L	1g	16,100g
XS32000L	TBA	Balance, XS32000L	1g	32,100g

Contact us ...



ADELAIDE: 08 8186 0523 rowesa@rowe.com.au  
 BRISBANE: 07 3376 9411 roweqld@rowe.com.au  
 MELBOURNE: 03 8795 7771 rowevic@rowe.com.au  
 PERTH: 08 9302 1911 rowewa@rowe.com.au

[www.rowe.com.au](http://www.rowe.com.au)

## Balances - Precision

Item	Rowe Code	Description	Readability	Weighing Capacity
11139590	IB3508	PB153-S/FACT Precision Balance with FACT	0.001g	151g
11139594	IB3533	PB303-S/FACT Precision Balance with FACT	0.001g	310g
11139598	TBA	PB303-SDR/FACT Prec Bal with FACT DeltaRange®	0.001g / 0.01g	60 / 310 g
11139602	IB1925	PB403-S/FACT Precision Balance with FACT	0.001g	410g
11139606	IB3733	PB503-S/FACT Precision Balance with FACT	0.001g	510g
11139610	TBA	PB602-S/FACT Precision Balance with FACT	0.01g	610g
11139614	IB3742	PB1502-S/FACT Precision Balance with FACT	0.01g	1510g
11139618	IB3743	PB3002-S/FACT Precision Balance with FACT	0.01g	3100g
11139622	TBA	PB3002-SDR/FACT Prec Bal with FACT DeltaRange®	0.01g / 0.1g	600 / 3100 g
11139626	IB1979	PB4002-S/FACT Precision Balance with FACT	0.01g	4100g
11139630	TBA	PB1501-S/FACT Precision Balance with FACT	0.1g	1510g
11139634	TBA	PB3001-S/FACT Precision Balance with FACT	0.1g	3100g
11139638	IB3598	PB5001-S/FACT Precision Balance with FACT	0.1g	5100g
11139642	IB1958	PB8001-S/FACT Precision Balance with FACT	0.1g	8100g
11139646		PB8000-S/FACT Precision Balance with FACT	1g	8100g



Classic Plus precision balances PB-S/FACT



Classic precision balances AB



AB265-S	IB3546	Balance, AB265-S	0.01mg/ 0.1mg	61g/220g
AB54-S	IB3524	Balance, AB54-S	0.1mg	51g
AB104-S	IB3512	Balance, AB104-S	0.1mg	110g
AB204-S	IB3527	Balance, AB204-S	0.1mg	220g
AB304-S	IB3567	Balance, AB304-S	0.1mg	320g

AB-S includes RS 232 interface and internal calibration

Contact us ...



ADELAIDE: 08 8186 0523 rowesa@rowe.com.au  
 BRISBANE: 07 3376 9411 roweqld@rowe.com.au  
 MELBOURNE: 03 8795 7771 rowevic@rowe.com.au  
 PERTH: 08 9302 1911 rowewa@rowe.com.au

[www.rowe.com.au](http://www.rowe.com.au)



### Classic precision balances PB-L and SB



Item	Rowe Code	Description	Readability	Weighing Capacity
11139526	IB3506	PB153-L Precision Balance	0.001g	151g
11139530	IB3544	PB303-L Precision Balance	0.001g	310g
11139534	TBA	PB303-LDR Precision Balance Delta Range	0.001g / 0.01g	60 / 310 g
11139538	IB1686	PB602-L Precision Balance	0.01g	610g
11139542	TBA	PB1502-L Precision Balance	0.01g	1510g
11139546	IB3736	PB3002-L Precision Balance	0.01g	3100g
11139550	TBA	PB3002-LDR Precision Balance Delta Range	0.01g / 0.1g	600 / 3100 g
11139554	TBA	PB1501-L Precision Balance	0.1g	1510g
11139558	TBA	PB3001-L Precision Balance	0.1g	3100g
11139562	IB3596	PB5001-L Precision Balance	0.1g	5100g
11139566	IB1955	PB8001-L Precision Balance	0.1g	8100g
11139570	IB3528	PB8000-L Precision Balance	1g	8100g
SB8001	IB1953	Balance, SB8001	0.1g	8100g
SB12001	TBA	Balance, SB12001	0.1g	12100g
SB16001	IB2200	Balance, SB16001	0.1g	16100g
SB16001DR	IB3746	Balance, SB16001 Delta Range	0.1g / 1g	3200/ 16100g
SB24001DR	TBA	Balance, SB24001 Delta Range	0.1g / 1g	4800/ 24100g
SB32001DR	IB2036	Balance, SB32001 Delta Range	0.1g / 1g	6400/ 32100g
SB8000	TBA	Balance, SB8000	1g	8100g
SB16000	IB2210	Balance, SB16000	1g	16100g
SB32000	IB2212	Balance, SB32000	1g	32100g

### Classic Light Precision Balances PL



Optional RS232 interface and Rechargeable battery. Calibration requires external weight (not included).

PL83-S/00	TBA	Balance	0.001 g	81 g
PL83-S/01	TBA	Balance with RS232	0.001 g	81 g
PL303/00	TBA	Balance	0.001 g	310 g
PL303/01	TBA	Balance with RS232	0.001 g	310 g
PL202-S/00	IB3543	Balance	0.01g	210 g
PL202-S/01	TBA	Balance with RS232	0.01g	210 g
PL602-S/00	IB1687	Balance	0.01g	610 g
PL602-S/01	TBA	Balance with RS232	0.01g	610 g
PL1502-S/00	TBA	Balance	0.01g	1510 g
PL1502-S/01	TBA	Balance with RS232	0.01g	1510 g
PL3002/01	TBA	Balance with RS232	0.01 g	3100 g
PL601-S/00	IB1680	Balance	0.1g	610 g
PL601-S/01	IB3716	Balance with RS232	0.1g	610 g
PL1501-S/00	IB2099	Balance	0.1g	1510 g
PL1501-S/01	TBA	Balance with RS232	0.1g	1510 g
PL3001-S/00	IB1920	Balance	0.1g	3100 g
PL3001-S/01	TBA	Balance with RS232	0.1g	3100 g
PL6001-S/00	TBA	Balance	0.1g	6100 g
PL6001-S/01	TBA	Balance with RS232	0.1g	6100 g
PL6000-S/00	IB3741	Balance	1 g	6100 g
PL6000-S/01		Balance with RS232	1 g	6100 g

Contact us ...



ADELAIDE: 08 8186 0523 rowesa@rowe.com.au  
 BRISBANE: 07 3376 9411 roweqld@rowe.com.au  
 MELBOURNE: 03 8795 7771 rowevic@rowe.com.au  
 PERTH: 08 9302 1911 rowewa@rowe.com.au

[www.rowe.com.au](http://www.rowe.com.au)

## Balances - Precision



### Explorer Pro

Ohaus® Explorer Pro Precision Top Loading balances offer a full range of capacities and functions to meet your most demanding requirements.

Item	Rowe Code	Readability	Weighing Capacity
EP413D	TBA	0.001 g / 0.01 g	100 g / 410 g
EP413DC	TBA	0.001 g / 0.01 g	100 g / 410 g
EP213	TBA	0.001 g	210 g
EP213C	TBA	0.001 g	210 g
EP413	IB7040	0.001 g	410 g
EP413C	IB4010	0.001 g	410 g
EP612	TBA	0.01 g	610 g
EP612C	TBA	0.01 g	610 g
EP613	TBA	0.001 g	610 g
EP613C	IB7461	0.001 g	610 g
EP4102D	TBA	0.01 g / 0.1 g	1000 g / 4100 g
EP4102DC	TBA	0.01 g / 0.1 g	1000 g / 4100 g
EP2102	IB7405	0.01 g	2100 g
EP2102C	IB7420	0.01 g	2100 g
EP4101	TBA	0.1 g	4100 g
EP4101C	TBA	0.1 g	4100 g
EP4102	IB7410	0.01 g	4100 g
EP4102C	IB7425	0.01 g	4100 g
EP6101	IB7436	0.1 g	6100 g
EP6101C	TBA	0.1 g	6100 g
EP6102	IB7433	0.01 g	6100 g
EP6102C	IB7434	0.01 g	6100 g
EP8101	TBA	0.1 g	8100 g
EP8101C	TBA	0.1 g	8100 g
EP12001	TBA	0.1 g	12000 g
EP12001C	TBA	0.1 g	12000 g
EP22001	TBA	0.1 g	22000 g
EP22001C	TBA	0.1 g	22000 g
EP32001	TBA	0.1 g	32000 g
EP32001C	TBA	0.1 g	32000 g



Contact us ...



ADELAIDE: 08 8186 0523 rowesa@rowe.com.au  
 BRISBANE: 07 3376 9411 roweqld@rowe.com.au  
 MELBOURNE: 03 8795 7771 rowevic@rowe.com.au  
 PERTH: 08 9302 1911 rowewa@rowe.com.au

[www.rowe.com.au](http://www.rowe.com.au)

## Balances - Precision



### Adventurer SL

Ohaus® Adventurer Precision Top Loading balances offer a full range of capacities and functions to meet your most demanding requirements.



Item	Rowe Code	Readability	Weighing Capacity
N10330	TBA	0.002 g	32 g
N12120	TBA	0.01 g	210 g
N14120	TBA	0.01 g	410 g
N06112	TBA	0.1 g	610 g
N18110	TBA	0.1 g	810 g
N1B110	IB7225	0.1 g	2100 g
N1D110	TBA	0.1 g	4100 g
N0F1A2	TBA	1 g	6100 g
N1H110	TBA	0.5 g	8100 g
N1H110	TBA	0.5 g	8100 g

### Adventurer



AR1530	IB7110	0.001 g	150 g
AR3130	IB7130	0.001 g	310 g
AR5120	TBA	0.01 g	510 g
ARA520	IB7124	0.01 g	1500 g
ARC120	IB3778	0.01 g	3100 g
ARD110	TBA	0.1 g	4100 g

### Navigator

The Ohaus® Navigator balance is a full-featured, compact, affordable, portable balance. The Navigator balance offers you capacities from 32g to 8100g.

(Note: Balances beginning with item # N1xxxx have Internal calibration and those with N0xxxx do not.)



N10330	TBA	0.002 g	32 g
N12120	TBA	0.01 g	210 g
N14120	TBA	0.01 g	410 g
N06112	TBA	0.1 g	610 g
N18110	TBA	0.1 g	810 g
N1B110	IB7225	0.1 g	2100 g
N1D110	TBA	0.1 g	4100 g
N0F1A2	TBA	1 g	6100 g
N1H110	TBA	0.5 g	8100 g

Contact us ...



ADELAIDE: 08 8186 0523 rowesa@rowe.com.au  
 BRISBANE: 07 3376 9411 roweql@rowe.com.au  
 MELBOURNE: 03 8795 7771 rowevic@rowe.com.au  
 PERTH: 08 9302 1911 rowewa@rowe.com.au

[www.rowe.com.au](http://www.rowe.com.au)

## Balances - Precision



The Voyager® Series was designed with the user in mind. Voyager Pro features an ultra-simple and innovative user interface. Key features include scroll down menus for simple menu navigation and on screen text prompts which guide you through Voyager Pro's application modes. Voyager® Pro's state of the art high-resolution display is easy to view and displays only the information you want to see. Plus, your investment is protected! Voyager® Pro is designed with a spill channel and includes an in-use cover. The Voyager® Pro is manufactured to ISO 9001:2000 quality assurance specifications.

Item	Rowe Code	Readability	Weighing Capacity
VP413DCN	TBA	0.001 g / 0.01 g	100 g / 410 g
VP213CN	TBA	0.001 g	210 g
VP413CN	TBA	0.001 g	410 g
VP612CN	TBA	0.01 g	610 g
VP4102DCN	TBA	0.01 g / 0.1 g	1000 g / 4100 g
VP2102CN	TBA	0.01 g	2100 g
VP4101CN	IB7600	0.1 g	4100 g
VP4102CN	TBA	0.01 g	4100 g
VP6101CN	TBA	0.1 g	6100 g
VP8101CN	TBA	0.1 g	8100 g

### Voyager Pro



AV53	IB7045	0.001 g	51 g
AV212	TBA	0.01 g	210 g
AV212C	TBA	0.01 g	210 g
AV213	TBA	0.001 g	210 g
AV213C	TBA	0.001 g	210 g
AV412	IB7182	0.01 g	410 g
AV412C	TBA	0.01 g	410 g
AV812	IB7185	0.01 g	810 g
AV812C	IB7183	0.01 g	810 g
AV2101	TBA	0.1 g	2100 g
AV2101C	TBA	0.1 g	2100 g
AV2102	IB7127	0.01 g	2100 g
AV2102C	IB7126	0.01 g	2100 g
AV4101	IB7190	0.1 g	4100 g
AV4101C	TBA	0.1 g	4100 g
AV8101	IB7580	0.1 g	8100 g
AV8101C	TBA	0.1 g	8100 g

### Adventurer Pro

The Ohaus® Adventurer™ Pro family combines unsurpassed value with solid reliability. Adventurer Pro balances are affordable, accurate and easy to use for all of your basic weighing needs. The Adventurer Pro family offers capacities ranging from 51g to 8100g.



Contact us ...



ADELAIDE: 08 8186 0523 rowesa@rowe.com.au  
 BRISBANE: 07 3376 9411 roweqld@rowe.com.au  
 MELBOURNE: 03 8795 7771 rowevic@rowe.com.au  
 PERTH: 08 9302 1911 rowewa@rowe.com.au

[www.rowe.com.au](http://www.rowe.com.au)

## Balances - Trade Approved



### Certificate of Approval No 6/4C/228

Characteristics of trade approved balances:

- All models feature internal calibration.
- The balance display shows the entire readability.
- Figures after the verification scale are shown in brackets.
- Internal calibration is automatically initiated on instrument switch on and 2 hours after switch on and every 2°C change and 150 hours.
- Only kilograms and grams are available units (except for -C jewellery models which have carats).
- External calibration is not available.
- Functions such as %, counting & dynamic weighing are fitted, but not approved for trade use.
- Instruments will carry appropriate labelling including certification mark.
- Instruments are sealed to prevent unauthorised calibration and tampering.
- All balances approved for use over temperature range 10-30°C except PB8000-S/A Class III version 10-40°C.

### Analytical Balances AB

Item	Rowe Code	Description	Weighing Capacity (Max) g	Verification "Scale Interval" (e) g	"Minimum capacity" "for trade use" (Min) g	Readability (d) g	Accuracy class
AB135-S/A	TBA	Balance, AB135-S/A	31 / 120	0.001 / 0.001	0.001	0.00001 / 0.0001	I
AB265-S/A	TBA	Balance, AB265-S/A	61 / 220	0.001 / 0.001	0.001	0.00001 / 0.0001	I
AB54-S/A	TBA	Balance, AB54-S/A	51	0.001	0.01	0.0001	I
AB104-S/A	TBA	Balance, AB104-S/A	110	0.001	0.01	0.0001	I
AB204-S/A	TBA	Balance, AB204-S/A	220	0.001	0.01	0.0001	I
AB304-S/A		Balance, AB304-S/A	320	0.001	0.01	0.0001	I



JB8001



JB3002



JB803

Contact us ...



ADELAIDE: 08 8186 0523 rowesa@rowe.com.au  
 BRISBANE: 07 3376 9411 roweqld@rowe.com.au  
 MELBOURNE: 03 8795 7771 rowevic@rowe.com.au  
 PERTH: 08 9302 1911 rowewa@rowe.com.au

[www.rowe.com.au](http://www.rowe.com.au)

## Balances - Trade Approved



### Precision Balances PB-S/SB

Item	"Rowe Code"	Description	Weighing Capacity (Max) g	Verification "Scale Interval" (e) g	"Minimum capacity" "for trade use" (Min) g	Readability (d) g	Accuracy class
PB153-S/A	TBA	Balance, PB153-S/A	151	0.01	0.02	0.001	II
PB303-S/A	TBA	Balance, PB303-S/A	310	0.01	0.02	0.001	II
PB303-SDR/A	TBA	Balance, PB303-S/A Delta	60 / 310	0.01 / 0.01	0.02	0.001 / 0.01	II
PB403-S/A	TBA	Balance, PB403-S/A	410	0.01	0.02	0.001	II
PB602-S/A	TBA	Balance, PB602-S/A	610	0.1	0.5	0.01	II
PB1502-S/A	TBA	Balance, PB1502-S/A	1510	0.1	0.5	0.01	II
PB3002-S/A	TBA	Balance, PB3002-S/A	3100	0.1	0.5	0.01	II
PB3002-SDR/A	TBA	Balance, PB3002-S/A Delta	600 / 3100	0.1 / 0.1	0.5	0.01 / 0.1	II
PB4002-S/A	TBA	Balance, PB4002-S/A	4100	0.1	0.5	0.01	II
PB1501-S/A	TBA	Balance, PB1501-S/A	1510	0.1	5	0.1	II
PB3001-S/A	TBA	Balance, PB3001-S/A	3100	0.1	5	0.1	II
PB5001-S/A	TBA	Balance, PB5001-S/A	5100	1	5	0.1	II
PB8001-S/A	TBA	Balance, PB8001-S/A	8100	1	5	0.1	II
PB8000-S/A	TBA	Balance, PB8000-S/A	8100	1	50	1	II
PB8000-S/A	TBA	Balance, PB8000-S/A	8100	1	20	1	III
JB3002-G/A	TBA	Balance, JB3002-G/A	3100	0.1	0.5	0.01	II
JB8001-G/A	TBA	Balance, JB8001-G/A	8100	1	5	0.1	II
JB203-C/A	TBA	Balance, JB203-C/A	51 255ct	0.01 0.01ct	0.02 0.02ct	0.001 0.001ct	II
JB803-C/A	TBA	Balance, JB803-C/A	162 810ct	0.01 0.01ct	0.02 0.02ct	0.001 0.001ct	II
JB1603-C/A	TBA	Balance, JB1603-C/A	320 1600ct	0.01 0.01ct	0.02 0.1ct	0.001 0.001ct	II I

### Jewellery Balances

JB3002-G/A	TBA	Balance, JB3002-G/A	3100	0.1	0.5	0.01	II
JB8001-G/A	TBA	Balance, JB8001-G/A	8100	1	5	0.1	II
JB203-C/A	TBA	Balance, JB203-C/A	51 255ct	0.01 0.01ct	0.02 0.02ct	0.001 0.001ct	II II
JB803-C/A	TBA	Balance, JB803-C/A	162 810ct	0.01 0.01ct	0.02 0.02ct	0.001 0.001ct	II II
JB1603-C/A	TBA	Balance, JB1603-C/A	320 1600ct	0.01 0.01ct	0.02 0.1ct	0.001 0.001ct	II I

Contact us ...



ADELAIDE: 08 8186 0523 rowesa@rowe.com.au  
 BRISBANE: 07 3376 9411 roweqld@rowe.com.au  
 MELBOURNE: 03 8795 7771 rowevic@rowe.com.au  
 PERTH: 08 9302 1911 rowewa@rowe.com.au

[www.rowe.com.au](http://www.rowe.com.au)

## Balances - Accessories

### Printers and Printer Accessories



Item	Rowe Code	Description
RS-P42	IT4220	Printer
BT-P42	TBA	Printer, Blue tooth
LC-P43	IB2025	Printer, simple dot matrix, LC interface for PR, SR, PG, SG, AB, PB, SB, HR73 & HG52
LC-P45	IB2030	Printer, versatile dot matrix, LC/RS232 interface for all MT Instruments
LC-PVolume	TBA	Pipette calibration printer for professional and standard balances
00072456	TBA	Paper, set of 5 rolls, for LC-P45, LC-P43, SQC14, GA42, RS-P42
00071328	TBA	Printer Paper, set of 4 rolls, for GA44
00065975	IB2020	Ink ribbon, 1 unit, for LC-P45, LC-P43, SQC14, GA42, RS-P42 (Must buy in Qty of 2)
11600388	TBA	Paper, adhesive, 1 roll, for LC-P45, LC-P43, SQC14, GA42, RS-P42
Anti-Theft		
590101	IB1999	Anti theft device, cable lock for AB-S/PB-S/SB/PG-S/SG/PR/SR/PL/AL
229175	TBA	Anti theft device, metal rod lock, pin with lock for AB/PB/PG/PR
11600361	IB1998	Anti theft device, steel cord, for XS/XP/AX
In-Use Covers		
11103681	IB2082	In use cover, for AB-S, PB-S
11100830	IB2073	In use cover, for front terminal, AX, MX, UMX
12102970	IB3230	In use cover, AL/PL
12102980	IB2014	In use cover, PL-S
11106870	TBA	In use cover, XS analytical & precision terminal
11100089	TBA	Dust cover for AX

### Interface Output Options

11103678	TBA	Interface RS232 for AB-S, PB-S, standard on models shipped after March 9, 2000
11100070	TBA	Interface RS232 for UMX, MX, AX balances
11100071	TBA	Interface LocalCan for UMX, MX, AX balances
11100072	TBA	Interface Mini Mettler for UMX, MX, AX balances
11132500	TBA	Interface, RS232 for XS (additional)
11132505	TBA	Interface, LC option, LocalCAN for XS
11132510	TBA	Interface Mini Mettler for XS & XP
11132520	TBA	Interface PS/2 for XS & XP
11132535	TBA	Interface Bluetooth for XS & XP
11132515	TBA	Interface Ethernet for XS & XP

### Software

237010	TBA	BalanceLink Software 3.5" Disk (Windows and NT)
--------	-----	---

Contact us ...



ADELAIDE: 08 8186 0523 rowesa@rowe.com.au  
 BRISBANE: 07 3376 9411 roweql@rowe.com.au  
 MELBOURNE: 03 8795 7771 rowevic@rowe.com.au  
 PERTH: 08 9302 1911 rowewa@rowe.com.au

[www.rowe.com.au](http://www.rowe.com.au)