

## XRF

### XRF Binder

For binding powder samples. 10g of sample is mixed with 2mL of XRF binder.

Rowe Code	Std. Pack Size
XB1000	2L



### XRF Cups

These cups allow the analysis of solutions by XRF, and fit all common makes of XRF instrument, including X-Unique II, PW2400, PW2404 and Axios, all Panalytical instruments.

Rowe Code	Description
XC1005	CUP, XRF, LIQUID (CLEAR)
XC1006	CUP, XRF, LIQUID (GREY)
XC1007	CUP, XRF, LIQUID (ORANGE)



### XRF Disc Tray

**Rowe Code: XX1050**

PVC Disc Tray 200mm (L) 45 (W) holds up to 42mm diameter discs.

### XRF Pellet Cups

Somar-caps are compressible thin-walled aluminium cups for briquetting powder samples. Simply fill the somar-cap with the powdered sample and briquette in a standard mould assembly - sample weighing is not required. The resultant pellet is reinforced with aluminium on all sides except the surface to be investigated. Somar-caps are chemically treated to eliminate adhesive to the mould assembly with provision on the backs for recording sample identification numbers.



Rowe Code	Diameter mm	Briquette Size mm	Std. Pack Size
XP1050	33	35	1500
XP1100	38	40	1000

### Other XRF related products:

Also see in this catalogue:  
Grinding Tablets  
Mylar film

### XRF/XRD Consumables

Rowe Scientific Pty Ltd supply XRF/XRD Consumables from the Chemplex, Somar and Spex ranges, and can supply a separate catalogue upon request.