



ADELAIDE
08 8186 0523
rowesa@rowe.com.au

BRISBANE
07 3376 9411
roweqld@rowe.com.au

MELBOURNE
03 8795 7771
rowevic@rowe.com.au

PERTH
08 9302 1911
rowewa@rowe.com.au

SYDNEY
02 9482 8433
rowensw@rowe.com.au

www.rowe.com.au

The **NEW** dispenser generation

The Labmax “airless” has an automatic air purging system. This avoids wrong values when dosing gas emitting chemicals, the lighter gas escapes automatically through the special ejection channel.

Contamination of bottle contents does not occur due to air purging being managed through the ejection cannula. The Labmax “airless” is made of high resistant materials and also offers a “HF” version for hydrofluoric acid dosing.

All dispenser-versions are fully autoclavable conformity certified. Each device includes a test certificate with corresponding identification number.

- Automatic air purging
- Safety-Lock
- Autoclavable
- Long-life cycle
- Adapters for bottles included
- HF version available

Both the “airless universal” and the “airless HF” dispensers are made of high-grade materials.

PART	AIRLESS UNIVERSAL	AIRLESS HF
Ejection Cannula	FEP	
Seal	PTFE	
Valve head	ETFE	
Cylinder	PTFE	
Plunger	Borosilicate glass	PP
Valve spring	Platinum Iridium	
Valve ball	Aluminium Oxide	Ruby
Valve seat	Borosilicate glass	Sapphire



Rowe Code	Volume Range	Description
DD1177	1-5mL,	AIRLESS DISPENSER
DD1176	10-50mL	SOLVENT, AIRLESS DISPENSER
DD1178	10-50mL,	AIRLESS DISPENSER
DD1175	2-10mL	SOLVENT, AIRLESS DISPENSER
DD1171	2-10mL,	AIRLESS DISPENSER
DD1179	25-100mL,	AIRLESS DISPENSER
DD1174	5-25mL	HF, AIRLESS DISPENSER
DL1106	5-25mL	SOLVENT, AIRLESS DISPENSER
DD1173	5-25mL,	AIRLESS DISPENSER