



# Rowe Scientific PTY LTD

*For accuracy and professionalism*



**ADELAIDE**  
08 8186 0523

**BRISBANE**  
07 3376 9411

**MELBOURNE**  
03 8795 7771

**PERTH**  
08 9302 1911

[www.rowe.com.au](http://www.rowe.com.au)

## XRF CUPS

These cups allow the analysis of solutions by XRF, and fit all common makes of XRF instruments, including X-Unique II, PW2400, PW2404, Axios, and many Panalytical instruments.

- Free trial samples available
- Very cost competitive.
- Avoids cross contamination between samples, as cups are disposable.
- Made from polypropylene, so chemically inert.
- Packaged under clean room conditions, so free from silica and other air borne particulates.



Rowe Code: XC1005    **Clear**    LID outside diameter: 45.7mm. BODY Inside diameter: 35.5mm  
Rowe Code: XC1006    **Grey**    LID outside diameter: 45.7mm. BODY Inside diameter: 35.5mm  
Rowe Code: XC1007    **Orange** LID outside diameter: 45.7mm. BODY Inside diameter: 35.5mm

To find out more or to acquire your **FREE** samples call your local **Rowe Scientific Pty. Ltd.** Office.

### Other XRF related products:

Grinding Tablets

Mylar film

XRF Binder for binding powder samples. 10g of sample is mixed with 2mL of XRF binder.

XRF Disc Tray PVC Disc Tray 200mm (L) 45 (W) holds up to 42mm diameter discs.

Rowe Scientific Pty Ltd supply XRF/XRD Consumables from the Chemplex, Somar and Spex ranges, and can supply a separate catalogue upon request.