



# Rowe Scientific PTY LTD

For accuracy and professionalism



## PRODUCT INFORMATION

ADELAIDE 08 8186 0523 - BRISBANE 07 3376 9411 - MELBOURNE 03 8795 7771 - PERTH 08 9302 1911 - www.rowe.com.au



# HQ40d multi-parameter digital meter

ROWE CODE IM6700



Plug and Play with IntelliCAL™ probes -pH, conductivity, LDO.

- Dual display of two different parameters/probes!
- A single meter measures pH, Conductivity, and LDO®
- Simultaneously measure two different parameters at once
- View information from two probes on one screen (up to two parameters)
- Transfer GLP/ISO data effortlessly to PC or flash stick/transfer entire datalog, or individual readings
- Includes Field Kit, which contains a protective glove, 2 probe holders, and sample cups in one organized portable kit
- Quick and easy swapping of parameter/probe without re-calibration
- Password-protected data for accurate, tamper proof reporting
- 500-point data log stores measurements, calibrations and check standards, offering complete traceability.

Comes with HQ40d Meter, manual, 4 AA batteries, power supply, USB/DC adapter and field kit (field kit includes protective meter glove, wrist strap, neck strap, 2 standard probe holders and 5 sample cups)



ROWE CODE IP5302



ROWE CODE IP5315

### NOTE: PROBES SOLD SEPARATELY.

*IntelliCAL™ Rugged Probes can reliably test in harsh environmental conditions. Move probes between meters without the need to re-calibrate. Super tough, stainless steel body and polymer sensor shroud. Choice of cable lengths- 5, 10, 15, or 30 meters. Track performance with calibration history. Probe alerts when re-calibration is needed. Ultimate traceability-easily obtain:- Time and Date stamp - Operator and Sample ID - Calibration history - Parameter - Reading. Factory calibrated probes - be up and running right out of the box*

*pH Probes Ranges pH 2 - 14 (Refillable PHC30101 has range pH 0 - 14)*

*To find out more on the probes available contact your Rowe Scientific Pty Ltd representative.*

### FEATURES

- HQ40d Dual Input
- AC/DC power or 4 AA batteries
- Measure pH, Conductivity and LDO
- Backlit screen

### PROBES

- IntelliCal probes
- Stores calibration
- LDO Luminescent Dissolved Oxygen
- Powder sachets / pillows

### BENEFITS

- Simultaneously measure using two different probes.
- Great for field use, 200+ hours battery life
- One metre to does all
- Great field feature in low light conditions

- World standards
- No calibration or re-calibration required
- No membranes to clean or replace
- No weighing or preparation of standards

### Related products to help you function

- Reagents
- Pipettes, Pipettors and Digital titrators
- Plasticware and general laboratory consumables.