



QUEENSLAND SPECIALS

100% AUSTRALIAN OWNED



Supplying scientific laboratories throughout Australia and overseas since 1987.

FREE TIMER

with orders placed before 31st May 2020.
Eligible for orders over \$150 dollars.
1 per client.



LABORATORY HYGIENE AND PERSONAL PROTECTION

Gloves

TouchNTuff® 92-600



CODE	DESCRIPTION	QTY	SIZE	PRICE
HG1891	Gloves disposable nitrile powder free Touch N Tuff	100	SML	\$25.50
HG1893	Gloves disposable nitrile powder free Touch N Tuff	100	MED	\$25.50
HG1895	Gloves disposable nitrile powder free Touch N Tuff	100	LGE	\$25.50
HG1897	Gloves disposable nitrile powder free Touch N Tuff	100	XLGE	\$25.50

Alcohol Wipes

GLR Wipes

Stocks arriving in June. Order Now!

CODE	DESCRIPTION	PRICE
HW0028	Wipes Alcohol "GLR", 42x14cm (70% Isopropyl Alcohol/30% Di Water)	\$9.95



Laboratory Coats

Long Sleeve Polypropylene

CODE	DESCRIPTION	QTY	PRICE
HC1775	Coat Laboratory Disposable White Polypropylene With No Pockets And Elastic Cuffs	50	\$104



Australian Made Surgical Face Masks

Surgical Mask

Surgical mask with ties or loops



ROWE CODE	DESCRIPTION	QTY	PRICE
HM0196	Mask, Surgical Face, Fluid Resistant With Ties, Green	6 x 50	\$240
HM0042	Mask, Surgical Face, Fluid Resistant With Loops	6 x 50	\$240

Manufactured to AS 4381:2002

Safety Glasses

Bollé protective spectacles



CODE	DESCRIPTION	PRICE
NG0050	Bollé Overlight II Safety Glasses Anti Scratch Anti Fog Clear Lens	\$10
NG0055	Bollé Prism Safety Glasses Anti Scratch Anti Fog Clear Lens	\$10
NG0063	Bollé Prism Seal Safety Glasses Anti Scratch Anti Fog Clear Lens	\$12
NG0037	Bollé Hustler Safety Glasses Anti Scratch Anti Fog Clear Lens Gloss Black Frame	\$23

Order by phone or online at www.rowe.com.au



FREE TIMER with orders placed before 31st May 2020.

Eligible for orders over \$150 dollars. 1 per client.

IT0240

Liquid Handling

DLAB LEVO Pipette Controller

- Precise pipetting control
- Robust and lightweight
- Compatible with most of the plastic and glass pipettes from 0.1 -100mL
- 3µm replaceable hydrophobic filter



CODE	DESCRIPTION	PRICE
DP0099	Pipette Controller LEVO 0.1-100mL Blue	\$67.00
DP0133	Pipette Controller LEVO 0.1-100mL Green	\$67.00
DP0132	Pipette Controller LEVO 0.1-100mL Red	\$67.00
DP0134	Pipette Controller LEVO 0.1-100mL Yellow	\$67.00

Bottle-top dispensers

Fortuna Optifix

- High quality materials and robust design for a long life-span
- Safe and efficient dispensing of dangerous chemicals



CODE	DESCRIPTION		PRICE
DL5300	Optifix Dispenser 2mL Basic	2mL	\$420.47
DL6000	Optifix Dispenser 2mL Safety	2mL	\$547.86
DL5350	Optifix Dispenser 5mL Basic	5mL	\$420.47
DL6030	Optifix Dispenser 5mL Safety	5mL	\$547.86
DL5160	Optifix Dispenser 5mL Solvent	5mL	\$420.47
DL5400	Optifix Dispenser 10mL Basic	10mL	\$420.47
DL5170	Optifix Dispenser 10mL Solvent	10mL	\$420.47
DL5450	Optifix Dispenser 30mL Basic	30mL	\$579.71
DL6060	Optifix Dispenser 30mL Safety	30mL	\$753.13
DL5500	Optifix Dispenser 50mL Basic	50mL	\$634.54
DL6100	Optifix Dispenser 50mL Safety	50mL	\$823.15
DL5180	Optifix Dispenser 50mL Solvent	50mL	\$634.54
DL6120	Optifix Dispenser 100mL Safety	100mL	\$1,168.06
DL5190	Optifix Dispenser 100mL Solvent	100mL	\$1,067.91

DLAB LEVO PLUS

The economical solution

CODE	DESCRIPTION	PRICE
DP0094	Pipette filler "LEVO PLUS", 0.1-100mL, With AC Adapter 17400032 And Filter 17000103	\$285

Aspiration Speed	25mL <7s (8 shift)
Dispensing Speed	Motor (8 shift)/Gravity
Battery	Lithium-ion
Battery Service Life	More than 8 Hours of intermittent Use
Charging Time	2-3 Hours
Pipette Types	Glass or Plastic Pipette (0.1-100mL), Pasteur Pipettes
Filter	0.45µm Hydrophobic
Weight	200g



- Large LCD display showing low battery warning and speed settings
- Select the desired coarse speed by pressing the button +/-

- Compatible with most plastic and glass pipettes from 0.1 -100mL
- Eight speeds for aspiration and dispensing different liquids
- Large LCD display showing low battery warning & speed settings
- Enables single-handed operation with minimum effort
- Light and ergonomic design providing easy usability
- High capacity Li-ion battery enables long runtime of operation
- Powerful pump fills a 25mL pipette in <7 seconds
- 0.45µm replaceable hydrophobic filter

Which BRAND™ dispenser do you need?

CODE	DESCRIPTION	Capacity	PRICE
DB0018	Dispensette S Analog-Adjustable 1-10mL x 0.2mL With Recirculation Valve	1-10mL	\$695
DD0036	Dispensette S Analog-Adjustable 10-100mL x 1.0mL With Recirculation Valve	10-100mL	\$1295
DD0017	Dispensette S Digital 1-10mL With Recirculation Valve	1-10mL	\$795
DD0024	Dispensette S Digital 5-50mL With Recirculation Valve	5-50mL	\$990





A 100% Australian owned company since 1987.

ROWE SCIENTIFIC PTY LTD

QUEENSLAND SPECIALS

Glassware

Laboratory Bottles

Borosilicate glass from 50mL to 10L



CODE	DESCRIPTION	Vol mL	PRICE
GB7855	Laboratory Bottle With Blue Screw Cap Borosilicate Glass 50mL	50	\$5.65
GB7910	Laboratory Bottle With Blue Screw Cap Borosilicate Glass 100mL	100	\$4.22
GB7955	Laboratory Bottle With Blue Screw Cap Borosilicate Glass 250mL	250	\$4.84
GB8005	Laboratory Bottle With Blue Screw Cap Borosilicate Glass 500mL	500	\$5.98
GB8055	Laboratory Bottle With Blue Screw Cap Borosilicate Glass 1L	1000	\$9.84
GB8105	Laboratory Bottle With Blue Screw Cap Borosilicate Glass 2L	2000	\$16.20
GB8155	Laboratory Bottle With Blue Screw Cap Borosilicate Glass 5L	5000	\$105.40
GB8195	Laboratory Bottle With Blue Screw Cap Borosilicate Glass 10L	10000	\$215.90

Cell Plasticware

SPL Culture Flasks

SPL culture flasks are ergonomically designed for easy handling, and minimise the chance of contamination. Short, wide and angled neck design enables easy access of serological pipettes and cell scrapers.



CODE	DESCRIPTION	PRICE
PC0329	T25 SPL Polystyrene with HPDE Filter Cap (vented) TC Treated Sterile. 40 Bags x 5 Flasks	\$129.59
PC0330	T75 SPL Polystyrene With HPDE Filter Cap (vented) TC Treated Sterile. 20 Bags x 5 Flasks	\$125.69
PC0331	T175 SPL Polystyrene With HPDE Filter Cap (vented) TC Treated Sterile. 8 Bags x 5 Flasks	\$100.55

Ask us about
Beakers and flasks

Ask us about
Submersible stirrers



Stirring



VELP MST MAGNETIC STIRRER

Small, cost efficient and effective magnetic stirrer for flasks up to 5 litres.

Rowe Code: IS0152

\$165



VELP ARE ALUMINIUM HOTPLATE STIRRER

Aluminium hotplate stirrer of flasks up to 15 litres and a maximum of 370°C heating capacity.

Rowe Code: IM0079

\$605



VELP AREX-6 DIGITAL PRO HOTPLATE STIRRER

Safe and powerful digital ceramic coated hot plate stirrer with timer. Complete with Pt100 temperature probe for direct temperature control of the sample.

Rowe Code: IM0557

\$1,365



Transferpette® S Starter Kits

CODE	DESCRIPTION	PRICE
DP0074	Pipettor Set Transferpette S, Midi Starter Kit, Consists Of: 3 Pipettes (20µL, 200µL & 1000µL Variable Type), 3 X Tipbox Filled, 3 Shelf/Rack Mounts	\$698.00
DP0075	Pipettor Set Transferpette S, Standard Starter Kit, Consists Of: 3 Pipettes (10µL, 100µL & 1000µL Variable Type), 3 X Tipbox Filled, 3 Shelf/Rack Mounts	\$714.96
DP0181	Pipettor Set Transferpette S, Pipette Package 1, Consists Of: 5 Pipettes (2,5µL, 10µL, 100µL, 200µL & 1000µL), 3 X Tipbox Filled, Bench-Top Rack For 6 Transferpette S	\$1,245.42



Offer Ends 30th June 2020 or while stocks last.



Order by phone or online at **www.rowe.com.au**
FREE TIMER with orders placed before 31st May 2020. Act quick!

Eligible for orders over \$150 dollars. 1 per client.

Water Analysis

HQ40d Portable pH, Conductivity, Dissolved Oxygen, ORP, and ISE Multi-Parameter Meter

Meter kit includes field case, 4 AA batteries, protective meter glove, power adapter, USB/DC power adapter for data transfer, quick-start guide, user manual, and documentation CD.

Rowe Code: IM6700



\$2,565

DR1900 Portable Spectrophotometer

DR 1900 Portable Spectrophotometer includes basic printed instrument manual, dust cover, 1 in. square glass sample cells, AA batteries and adapter set.

Rowe Code: IS0019



\$6,359

DR900 Handheld Colorimeter

DR 900 Colorimeter, two 1-inch glass sample cells marked at 10, 20 and 25 mL, two 1 cm plastic sample cells, 1 x 16mm COD/Test 'N Tube™ adapter, 4 AA alkaline batteries, manuals, USB cables.


Rowe Code: IC1581



\$3,325

MAINTENANCE - REPAIRS - CALIBRATION

SERVICING



SERVICING

Send your instruments to us and we will handle the rest.



CALIBRATION

For accuracy.



REPAIRS

To most broken instruments.



MAINTENANCE PLANS

Cost efficient service plans to keep your laboratory instruments operating.



Here are just some of the brands we can help service and repair:



www.rowe.com.au

South Australia & NT
 Ph: (08) 8186 0523
 rowesa@rowe.com.au

Queensland
 Ph: (07) 3376 9411
 roweqld@rowe.com.au

Victoria & Tasmania
 Ph: (03) 9701 7077
 rowevic@rowe.com.au

New South Wales
 Ph: (02) 9603 1205
 rowensw@rowe.com.au

Western Australia
 Ph: (08) 9302 1911
 rowewa@rowe.com.au