

# NDT 700 REFRIGERATOR AND FREEZER

**NULINE**  
REFRIGERATION

The NDT Spark Free contains a dual refrigeration system: a unit for refrigeration and another for freezer.

This locally manufactured cabinet offers quality 304 grade stainless steel interior and exterior as standard with the option of 316.

## FEATURES

- Top mounted system
- Heated cabinet facia to prevent condensation
- Auto defrost in 2°C Section
- Manual defrost in -20°C section
- Stainless steel interior / exterior
- Digital display with built in alarms
- Hi/Lo Readout




Refrigerator: 2°C to 8°C  
Freezer: -18°C to -22°C



Refrigerator: 386 L  
Freezer: 212L



870w x 750d x 2120h



# NDT 700 REFRIGERATOR AND FREEZER

**NULINE**  
REFRIGERATION

MODEL	REFRIGERATOR	FREEZER
Gross Capacity	386 L	212 L
External Dimensions in mm (w/d/h including castors or legs)	870 x 750 x 2120	
Internal Dimensions in mm (w/d/h)	770 x 725 x 885	685 x 535 x 580
Packed Dimensions in mm (w/d/h including pallet)	920 x 800 x 2250	
Gross Weight – kg (packed weight)	132	
Net Weight – kg	125	
Voltage	240V	
Temperature Range	2°C to 8°C	-18°C to -22°C
Ambient Temperature Range	10C° to 32C°	
Refrigerant	R134a	R404
Cooling System / Defrost	Static	Static
Exterior Cabinet Finish	Stainless Steel	Stainless Steel
Inner Liner Material	Stainless Steel	Stainless Steel
Door	Solid	Solid
Door Hinge	RH	RH
Door Handle	Lockable	Lockable
Type of Controller	Digital	Digital
Temperature Display	External	External
Data Logger	Optional	Optional
Door Ajar Alarms	Optional	Optional
High and Low Readout	Visual / Audible	Visual / Audible
Remote Alarm Contacts (BMS Connectivity for Alarms only) (BMS Connectivity for Temperatures)	Optional Optional	Optional Optional
Port Entry	Optional	Optional
Battery Back-Up for Controller	Optional	Optional
Power Requirement	10amp GPO	10amp GPO
Watts	960W	540W
Castors	Yes	
Adjustable Shelves	3	2
Usable Shelf Area in mm (w/d)	750 x 555	520 x 675



**South Australia & NT**  
Ph: (08) 8186 0523  
rowesa@rowe.com.au

**Queensland**  
Ph: (07) 3376 9411  
roweqld@rowe.com.au

**Victoria & Tasmania**  
Ph: (03) 9701 7077  
rowevic@rowe.com.au

**New South Wales**  
Ph: (02) 9603 1205  
rowensw@rowe.com.au

**Western Australia**  
Ph: (08) 9302 1911  
rowewa@rowe.com.au

[www.rowe.com.au](http://www.rowe.com.au)

