




**ROWE SCIENTIFIC** PTY LTD www.rowe.com.au

For accuracy and professionalism

Providing laboratory supplies to the scientific community across Australia since 1987.

We are proudly a 100% Australian owned company.



# Clean your lab glassware and apparatus with **QUANTUMCLEAN Concentrate**



Proudly manufactured in Australia by Rowe Scientific Pty Ltd



**Quantumclean.** For total cleaning (including radioactive decontamination) of all laboratory and industrial glassware, plastics, ceramics and stainless steel.

- **Ideal for use in ultrasonic cleaners**
- **Phosphate free**
- **Biodegradable**

**USE:** Supplied as a liquid concentrate for dilution with water (approx. 2-5% is adequate for most laboratory applications). Ideal for manual cleaning or for use in an ultrasonic cleaner.



Royal Flying Doctor Service

**PROUDLY SUPPORTING THE ROYAL FLYING DOCTOR SERVICE OF AUSTRALIA**

Your purchases of 1, 5 and 20L Quantumclean and 500mL Rowe pH buffer in the unique twin neck bottles will result in \$2 per item being sent to the RFDS to support this great Australian institution.

## [www.rowe.com.au](http://www.rowe.com.au)

**South Australia & NT**  
Ph: (08) 8186 0523  
rowesa@rowe.com.au

**Queensland**  
Ph: (07) 3376 9411  
roweqld@rowe.com.au

**Victoria & Tasmania**  
Ph: (03) 9701 7077  
rowevic@rowe.com.au

**New South Wales**  
Ph: (02) 9603 1205  
rowensw@rowe.com.au

**Western Australia**  
Ph: (08) 9302 1911  
rowewa@rowe.com.au