



**ADELAIDE**  
08 8186 0523  
rowesa@rowe.com.au

**BRISBANE**  
07 3376 9411  
roweqld@rowe.com.au

**HOBART**  
03 6272 0661  
rowetas@rowe.com.au

**MELBOURNE**  
03 8795 7771  
rowevic@rowe.com.au

**PERTH**  
08 9302 1911  
rowewa@rowe.com.au

**SYDNEY**  
02 9482 8433  
rowesw@rowe.com.au


**www.rowe.com.au**

## Volac Bottles

Polystop® **Dropping Bottles**, clear or amber glass, have Polystop® 3 piece stoppers (including a strong glass dropper, a durable vinyl teat and a polypropylene stopper). Full spares are available, packs of 10.

Polystop® **Dispensing Bottles** allow pipetting of measured volumes from a robust graduated glass dropping pipette.

Rowe Scientific Pty Ltd can now source all of your Volac dropping and dispensing bottle needs. To find out more please talk to your local Rowe Scientific consultant.



### Clear dispensing Bottles (1)

Rowe Code	Capacity	Dropper Volume
GB4608	100mL	1mL - 0.5mL
GB4630	250mL	8mL - 0.2mL
GB4609	500mL	10mL - 0.2mL

### Clear dropping Bottles (2)

Rowe Code	Capacity	Stopper (neck mm)
GB4627	50mL	16 / 16
GB4638	100mL	16 / 16

### Amber dispensing Bottles

Rowe Code	Capacity	Dropper Volume
GB4780	100mL	1mL - 0.5mL
GB4785	250mL	8mL - 0.2mL
GB4790	500mL	10mL - 0.2mL

### Amber dropping Bottles (3)

Rowe Code	Capacity	Stopper (neck mm)
GB4607	50mL	16 / 16
GB4639	100mL	16 / 16