

A 100% Australian owned company supplying scientific laboratories throughout Australia and overseas since 1987.



Some of the latest and most popular laboratory products available from Rowe Scientific Pty Ltd.



October 4 1957 Launch of Sputnik I, the first artificial satellite to orbit the Earth.



November 11 1930 Patent number US1781541 is awarded to Albert Einstein and Leó Szilárd for their invention, the Einstein refrigerator.



December 3 1967 First successful heart transplant carried out by Christiaan Barnard.

Source: en.wikipedia.org/wiki/Portal:History_of_science. Images Shutterstock.

Sterile Syringe Filter Membranes

Membranes available in a range of diameters and pore sizes. Membrane types:

PES (polyethersulfone), **PVDF** (polyvinylidene difluoride), **MCE** (mixed cellulose esters), **CA** (cellulose acetate), **CN** (cellulose nitrate), nylon, polypropylene **PTFE** (polytetrafluoroethylene), polycarbonate and glass-fibre prefilter.

Samples are available for trial. **FILTER-BIO**



ROWE CODE	DESCRIPTION	UOM	PRICE
FF0256	Syringe filter 13mm 0.22 micron, PTFE - hydrophobic gamma sterile, individually packed	100	\$119.70
FF0098	Syringe filter 25mm 0.22 micron PES, sterile, individually packed	50	\$54.60
FF0097	Syringe filter 25mm 0.45 micron PES, sterile, individually packed	50	\$54.60
FF0022	Syringe filter 25mm x 0.22 micron nylon sterile	50	\$87.15
FF0248	Syringe filter 30mm 0.22 micron polyvinylidene difluoride (PVDF) gamma sterile, individually packed	50	\$84.00
FS0049	Syringe filter 33mm 0.22 micron MCE acrylic housing, sterile	50	\$94.50
FS0046	Syringe filter 33mm 0.45 micron PVDF hydrophilic acrylic housing, sterile	50	\$94.50
FF0124	Syringe filter sterile 13mm 0.22 micron cellulose acetate	100	\$119.70
FF0128	Syringe filter sterile 25mm 0.45 micron cellulose acetate	50	\$63.00



- ✓ One stop shop
- ✓ 100% Australian
- ✓ Online ordering
- ✓ Phone ordering
- ✓ Account managers
- ✓ Leading brands
- ✓ Chemicals
- ✓ Consumables
- ✓ Equipment

Below are just a few leading brands we offer.



See our range of PPE products on pages 5 & 6. Note: These products are currently in high demand and limited supply.



TISSUE CULTURE

Culture flasks from SPL Life Sciences are ergonomically designed for easy handling, while achieving minimal contamination during cell culture. Short, wide and angled neck design enables easy access of serological pipettes and cell scrapers.

Flasks are divided into 6 types, depending on the culture volume and the cap type. Each flask can be provided with plug or filter caps.



CELL FLASK GROWTH AREA
T25 25cm² **T75** 75cm² **T175** 175cm²

	ROWE CODE	DESCRIPTION	UOM	PRICE
BIOREACTOR	PB0445	Bioreactor with 4 adjustable levels of air circulation, 15mL, sterile, 30 x 3pk	90	\$77.57
	PB0446	Bioreactor with 4 adjustable levels of air circulation, 50mL, sterile, 30 x 3pk	90	\$86.19
CELL FLASK	PC0329	Cell culture flask T25, SPL, 25cm ² , polystyrene with HDPE filter cap (vented), TC treated, sterile, 40 bags x 5 flasks, DNase/RNase free, non-cytotoxic, non-pyrogenic	200	\$113.25
	PC0330	Cell culture flask T75, SPL, 75cm ² , polystyrene with HDPE filter cap (vented), TC treated, sterile, 20 bags x 5 flasks, DNase/RNase free, non-cytotoxic, non-pyrogenic	100	\$109.85
	PC0331	Cell culture flask T175, SPL, 175cm ² , polystyrene with HDPE filter cap (vented), TC treated, sterile, 8 bags x 5 flasks	40	\$87.85
	PC0388	Cell culture flask T25, SPL, 25cm ² , plug cap, surface treated, polystyrene/HDPE sterile, 40 x 5pk	200	\$111.40
	PC0389	Cell culture flask T75, SPL, 75cm ² , plug cap, surface treated, polystyrene/HDPE sterile, 20 x 5pk	100	\$104.45
	PC0391	Cell culture flask T175, SPL, 175cm ² , plug cap, surface treated, polystyrene/HDPE sterile, 8 x 5pk	40	\$84.65
	PC0385	Cell culture flask T25, SPL, 25cm ² , polystyrene with HDPE filter cap (vented), non TC treated, sterile, 40 bags x 5 flasks, DNase/RNase free, non-cytotoxic, non-pyrogenic	200	\$116.50
	PC0390	Cell culture flask T75, SPL, 75cm ² , filter cap, untreated, polystyrene/HDPE sterile, 20 x 5pk	100	\$109.80
	PC0392	Cell culture flask T175, SPL, 175cm ² , filter cap, untreated (for suspension culture) polystyrene/HDPE sterile, 8 x 5pk	40	\$99.45
CELL STRAINER	PC0348	Cell strainers, 40 um, blue, SPL, sterile, PP/nylon fits 50mL conical tubes	50	\$63.05
	PC0349	Cell strainers, 70 um, natural, SPL, sterile, PP/nylon fits 50mL conical tubes	50	\$60.45
	PC0350	Cell strainers, 100um, yellow, SPL, sterile, PP/nylon fits 50mL conical tubes	50	\$64.40
CELL CULTURE	PC0332	Cell culture dish 35x10mm, SPL, TC treated, external grip, sterile, polystyrene, 50 bags of 10, DNase/RNase free, non-cytotoxic, non-pyrogenic	500	\$93.10
	PC0333	Cell culture dish 60x15mm, SPL, TC treated, external grip, sterile, polystyrene, 50 bags of 10, DNase/RNase free, non-cytotoxic, non-pyrogenic	500	\$112.75
	PC0334	Cell culture dish 90x20mm, SPL, TC treated, external grip, sterile, polystyrene, 20 bags of 10, DNase/RNase free, non-cytotoxic, non-pyrogenic	200	\$68.80
	PC0335	Cell culture dish 150x25mm, SPL, TC treated, sterile, polystyrene, 12 bags of 10, DNase/RNase free, non-cytotoxic, non-pyrogenic	120	\$128.15
	PC0337	Cell culture plate, 12 well, SPL, 85.4X127.6mm, flat bottom, TC treated, sterile, polystyrene, DNase/RNase free, non-cytotoxic, non-pyrogenic	50	\$62.35
	PC0381	Cell culture plate, 24 well, SPL, 85.4X127.6mm, flat bottom, TC treated, sterile, polystyrene, individually wrapped, DNase/RNase free, non-cytotoxic, non-pyrogenic	50	\$61.75
	PC0338	Cell culture plate, 96 well, SPL, 85.4X127.6mm, flat bottom, TC treated, sterile, polystyrene, DNase/RNase free, non-cytotoxic, non-pyrogenic	50	\$69.85
	PC0339	Cell culture plate, 96 well, SPL, 85.4X127.6mm, round bottom, TC treated, sterile, polystyrene, DNase/RNase free, non-cytotoxic, non-pyrogenic	50	\$69.85

	ROWE CODE	DESCRIPTION	UOM	PRICE
ERLENMEYER	PE0135	Erlenmeyer cell culture flask, 250mL, SPL, filter cap, sterile, PC/PP (body/cap), DNase/RNase free	8	\$48.70
	PE0136	Erlenmeyer cell culture flask, 500mL, SPL, filter cap, sterile, PC/PP (body/cap), DNase/RNase free	8	\$62.65
	PE0137	Erlenmeyer cell culture flask, 1000mL, SPL, filter cap, sterile, PC/PP (body/cap), DNase/RNase free	4	\$51.10
	PE0138	Erlenmeyer cell culture flask, 2000mL, SPL, filter cap, sterile, PC/PP (body/cap), DNase/RNase free	6	\$161.50
SCRAPERS	PC0342	Cell scraper, length 300mm, SPL, blade width 20mm, sterile, polystyrene/LDPE, DNase/RNase free, non-cytotoxic, non-pyrogenic	100	\$54.85
PCR TUBE	PP0128	PCR tube, 0.2mL, single, dome cap, SPL, PP, 2 bags of 500	1000	\$60.50
SEROLOGICAL PIPETTES	PS0332	Serological pipette 1mL, SPL, yellow, individual, sterile, 1/200/800, DNase/RNase free, non-cytotoxic, non-pyrogenic	800	\$129.45
	PS0333	Serological pipette 2mL, SPL, green, individual, sterile, 1/150/600, DNase/RNase free, non-cytotoxic, non-pyrogenic	600	\$120.25
	PS0334	Serological pipette 5mL, SPL, blue, individual, sterile, 1/100/400, DNase/RNase free, non-cytotoxic, non-pyrogenic	400	\$92.20
	PS0335	Serological pipette 10mL, SPL, orange individual, sterile, 1/100/400, DNase/RNase free, non-cytotoxic, non-pyrogenic	400	\$99.95
	PS0336	Serological pipette 25mL, SPL, red individual, sterile, 1/50/200, DNase/RNase free, non-cytotoxic, non-pyrogenic	200	\$76.00
	PS0337	Serological pipette 50mL, SPL, purple, individual, sterile, 1/40/160, DNase/RNase free, non-cytotoxic, non-pyrogenic	160	\$235.45
RODAC PLATE	PR0148	Rodac plate, 61x16.7mm, SPL, growth area 25cm ² , sterile PS, 25 bags of 20	500	\$87.85
SPREADER BAR	PS0338	Spreader bar, sterile PP, SPL, 50 bags of 10, DNase/RNase free	500	\$48.50

Serological pipettes, SPL, individual, sterile, DNase/RNase free, non-cytotoxic, non-pyrogenic



Erlenmeyer cell culture flasks, SPL, filter cap, sterile, PC/PP (body/cap), DNase/RNase free



PR0148 Rodac plate, 61x16.7mm, SPL, growth area 25cm², sterile PS, 25 bags of 20



PP0128 PCR tube, 0.2mL, single, dome cap, SPL, PP, 2 bags of 500



SPL Life Sciences provides customers with a wide range of high quality cell culture products, from basic cell culture vessels to more advanced platforms and accessories.

All SPL cell culture products meet international guidelines/regulations, such as ISO (International Organization for Standardization), ANSI (American National Standard Institute) and USP (United States Pharmacopeia, class VI). Polymer resins used to manufacture all products have followed international guidelines such as CONEG, SARA, Reach EC/1907/2006, RoHS 2002/95/EC. Biological tests (endotoxin, cytotoxicity, DNase/RNase, Human DNA) as well as cell attachment tests are regularly performed for quality control.





OVERSTOCKED



Be quick, don't miss the chance to save on selected BRAND liquid handling products.
ONLY WHILE STOCKS LAST.



ROWE CODE	DESCRIPTION	UOM	AVAILABLE	LIST	SPECIAL
DB0014	Dispensette S digital 2.5-25mL with recirculation valve (Analog version recommended for concentrated acids)	EA	3	1056.20	701.50
DD0025	Dispensette S Digital 5-50mL with recirculation valve (Analog version recommended for concentrated acids)	EA	2	1124.80	833.00
DL4981	Dispensette S organic analog-adjustable 0.5-5mL with recirculation valve	EA	3	748.91	595.30
DP1200	Pipettor Transferpette S 0.1-1uL adjustable volume	EA	6	357.12	240.30
DB0017	Dispensette S organic analog-adjustable 2.5-25mL with recirculation valve	EA	2	952.25	661.50
DP0197	Pipette Transferpette-8 15-300uL electronic 8 channel	EA	1	1778.58	1,317.00
DD0060	Dispensette S trace analysis 1-10mL TA valve spring with safetyprime recirculation valve (Not suitable for HF/NaOH)	EA	1	1666.84	1,096.90
GD0028	Dispensette S trace analysis 1-10mL Pt/Ir valve spring without safetyprime recirculation valve. Suitable for HF	EA	1	1565.87	1,053.80
DP0145	Pipettor Transferpette electronic 50-1000uL	EA	1	880.43	629.30
DD0105	Dispensette S organic, replacement discharge tube, 5 & 10mL with recirculation valve	EA	3	237.18	142.70



Transferpette® S
New and improved model, same German quality.

BRAND Transferpette®S adjustable pipettes CE-IVD, DE-M package.
Available for immediate delivery.

- Short stroke and low operating forces
- Universal tip cone
- Completely autoclavable at 121°C
- One-hand operation for right and left handed users
- Can be calibrated easily and without tools thanks to the Easy Calibration technology
- Enhanced slim shaft design for sampling from narrow vessels (centrifuge tubes)
- Operating temperature from +15°C to +40°C
- Vapour pressure up to 500 mbar
- Able to work with a liquid viscosity up to 260 mPa

Transferpette® S package includes:

5 x pipettes (0.1-2.5 µl, 0.5-10 µl, 2-20 µl, 20-200 µl, 100-1000 µl), 3 x TipBox pipette tips (0.5-20 µl, 2-200 µl, 50-1000 µl), benchtop rack which accommodates 6 pipettes.




DP0263
\$1295
Package Deal

HAIER range of Ultra Low Temperature freezers

State of the art units with the J model offering low power consumption, the BP series with extra low power consumption (via a variable speed compressor), and for your critical samples, the S series which have a back-up cooling system.





Because safety is paramount.

- ✓ Safety glasses
- ✓ Masks
- ✓ Gloves
- ✓ Face shields
- ✓ Hand sanitiser
- ✓ Alcohol wipes

PPE products are currently in high demand and limited supply, for this reason we ask you check with your local Rowe Scientific state office to discuss your supply needs.

If you do not see what you are seeking here, please feel free to contact a Rowe Scientific consultant about the full range of PPE items we have available.

Surgical face mask with ear loops

HM0210
50 per pack

Ultra Health surgical face mask with Anti-fog nasal strip and integrated earloops.

- Surgical Mask 3 ply with ear loops
- Disposable, latex free
- Level 2 (Sub Micron)
- 17.5x9.5cm Adult Masks
- Conforms to AS4381:2015
- ARTG LIST: 296802



KN95 face mask with ear loops Splash resistant Type IIR

HM0212
25 per pack

This disposable white KN95 face mask offers a comfortable user experience which is soft and easy to breathe through.

- Made from non-woven, static cotton, melt-blown fabric
- Individually packed
- Professionally tested in Australia
- 3 layers of protection
- 99.9% bacterial filtration efficiency
- TGA listed (ARTG #336093)



Hand Sanitiser W.H.O. Hand rub formulation 1

Rowe Scientific manufactured personal hand sanitiser in accordance with the World Health Organisation formulation 1.

Available in 3 different pack sizes for your convenience.

ROWE CODE	Pack Size
CH2261	16 x 250 mL
HH0030	500 mL
HH0032	10 L



Disposable Powder free Nitrile Gloves

200 per pack

Niclean 4.2 gives users exceptional levels of comfort and dexterity thanks to its softer, lightweight formulation while still providing the protections expected of nitrile gloves. Micro textured fingertips help to improve grip in wet or dry



Bollé clear over-glasses anti-fog face shields

HF0118

This Bollé face shield protects against splashes and droplets. Designed to fit the first responders' needs, it is simple to use, economical, protective, and can be worn over glasses.

- Anti-fog coating on both sides
- Adapted for use with prescription glasses
- Protection against splashes & fluids
- Elastic headband (one size fits all)
- XL size visor: 22 x 33cm
- Anti-splash tested to AS/NZS 1337.1



Alcohol wipes GLR

HW0028
75 wipes

GLR wipes are ideal for cleaning laboratory equipment, medical equipment, benchtops and office equipment. Containing 70% Isopropyl alcohol.





CRYOVIALS AND STORAGE



Cryovials are designed for storage of cells, specimens and solutions. Two different cap types, external & internal are available. Storable in ultra low temperature freezers.

- Use only in vapour-phase liquid nitrogen
- Distinctive external & internal cap design.
- External cap: External thread of the body fits perfectly into internal thread of the cap in helical form
- Internal cap: Internal thread of the body fits perfectly into the external thread of the cap in helical form
- Self-standing bottom
- Packed in a re-sealable zip-lock pack



Cryovial 1.8mL with external or internal cap thread.



Cryochill coders. Available in blue, green, red, white and yellow.

ROWE CODE	DESCRIPTION	UOM	PRICE
PC0344	Cryovial 1.8mL, external cap, without silicone washer, SPL, sterile PP, 10 bags of 50	500	\$119.70
PC0345	Cryovial 1.8mL, internal cap, with silicone washer, SPL, sterile PP, 10 bags of 50, DNase/RNase free, non-cytotoxic, non-pyrogenic	500	\$128.85
PC0346	Cryovial 5mL, external cap, with silicone washer, SPL, sterile PP, 10 bags of 50	500	\$199.80
PC0347	Cryovial 5mL, internal cap with silicone washer, SPL, sterile PP, 10 bags of 50	500	\$199.80
PV0017	Vial cryogenic (Cryochill) 1.8mL internal thread, 47mm H (with cap) 12mm OD screw cap with silicone o-ring, starfoot base sterile	500	\$139.80
PV0009	Vial cryogenic (Cryochill) 4.5mL internal threaded with screw cap silicon washer starfoot base sterile	10	\$5.50
PV1025	Vial cryogenic (Cryovial) 1.8mL 47mm high (with cap) 12mm OD 14mm cap OD, polypropylene with screw cap, sterile	500	\$98.50
PV1030	Vial cryogenic (Cryovial) 4.5mL 90mm H (with cap) 12mm OD polypropylene with screw cap, sterile	250	\$74.50

Storage Box for 2" Cryotubes

Also available in 64-place, 81-place, and 100-place sizes, with assorted colour options to organise and colour code your vials in your laboratory.



Freezer box, 2" cryotube storage, cardboard 81 place (9 x 9) green, 134 mm L x 134mm W x 51mm H

ROWE CODE	DESCRIPTION	UOM	PRICE
HF0085	Freezer box, 2" cryotube storage 100 place (10x10) with divider, 134 mm L x 134mm W x 51mm H, cardboard green	EA	\$1.85
HF0063	Freezer box, 2" cryotube storage, cardboard 81 place (9 x 9) orange, 134 mm L x 134mm W x 51mm H	EA	\$1.85

PIPETTE TIPS



Maxipense™ ultra low retention Pipette tips are available with standard and filter tip configuration.

Tarsons Maxipense™ ultra low retention surface is achieved by treating the polypropylene at the molecular level, modifying the surface of the polymer to have a very low coefficient of friction.

Pack size: 10 racks x 96 tips (all graduated)
Brand Compatibility: Includes Eppendorf®, Biohit® (Boeco), Finnpiptette®, Gilson®, Rainin®, Nishiryō®

Universal Tip Fit:
TARSONS tips are engineered to fit all popular pipettes.

Material:
Tarsons tips are manufactured from Medical Grade Resin. This ensures high quality products only.

Thin Walled Tips:
Discover thin wall tips, built with the latest injection moulding processes. The tips are high precision with excellent clarity as a result.

Highly Polished Moulds:
For extra smooth surfaces that reduce sample retention.

Rigorous Quality Control
All Tarsons tips are inspected by automated quality control procedures. The products are free of detectable DNase, RNase and pyrogen. Sterile tips are validated according to ISO 11137-2 regulations.

Manufacture:
Automated moulding and assembling are carried out in clean room conditions - this eliminates all chances of contamination.

Ultra Low Sample Retention
Improved flow dynamics ensures smooth flow of sticky liquids. Tarsons incorporate extremely hydrophobic surfaces for ultra low retention of samples.


ROWE CODE	DESCRIPTION	UOM	PRICE
PT0213	Tip pipette filter 10uL graduated racked sterile polypropylene, autoclavable	960	\$59.90
PT0214	Tip pipette filter 10uL XL / 20uL graduated racked sterile polypropylene, autoclavable	960	\$59.90
PT0218	Tip pipette filter 20uL graduated racked sterile polypropylene, autoclavable	960	\$59.90
PT0211	Tip pipette filter 100uL graduated racked sterile polypropylene, autoclavable	960	\$59.90
PT0216	Tip pipette filter 200uL graduated racked sterile polypropylene, autoclavable	960	\$59.90
PT0191	Tip pipette filter 1000uL (1mL) graduated racked sterile polypropylene (10 racks x 96 tips) to suit Eppendorf, Finnpiptette, Gilson, Nishiryō, Rainin	960	\$59.90
DT0025	Tip, pipette, 10uL, universal filter Maxipense ultra low retention graduated racked sterile (10 racks x 96 tips) to suit Biohit (Boeco), Eppendorf, Finnpiptette, Gilson, Nishiryō, Rainin	960	\$62.85
DT0026	Tip, pipette, 20uL, length 53mm, universal filter Maxipense ultra low retention graduated racked sterile (10 racks x 96 tips) to suit Biohit (Boeco), Eppendorf, Finnpiptette, Gilson, Nishiryō, Rainin	960	\$62.85
DT0028	Tip, pipette, 100uL, universal filter Maxipense ultra low retention graduated racked sterile (10 racks x 96 tips) to suit Biohit (Boeco), Eppendorf, Finnpiptette, Gilson, Nishiryō, Rainin	960	\$62.85
DT0029	Tip, pipette, 200uL, universal filter Maxipense ultra low retention graduated racked sterile (10 racks x 96 tips) to suit Biohit (Boeco), Eppendorf, Finnpiptette, Gilson, Nishiryō, Rainin	960	\$62.85
DT0031	Tip, pipette, 1000uL (1mL), universal filter Maxipense ultra low retention racked sterile (10 racks x 96 tips) to suit Eppendorf, Finnpiptette, Gilson, Nishiryō, Rainin	960	\$62.85

Why choose Tarsons tips?

TARSONS filter tips have a hydrophobic polyethylene filter which protects the sample from cross contamination that occurs due to aerosol transfer. Designed to be free of detectable DNase, RNase, and pyrogen.

The tips are a universal fit and are compatible with most pipettes on the market. They contain no additives which can interfere with your sample and the Maxipense™ filters are hydrophobic which allows maximum recovery of your samples.

Ideal for use in Molecular Biology, Microbiology, Cell Culture and Radioactive industries or similar applications.



DETERGENTS



CQ1050 CD1240 CD1430

ROWE CODE	DESCRIPTION	UOM	PRICE
CD1431	Detergent pyroneg liquid 6 x 2L bottles	12L	\$269.90
CD1430	Detergent pyroneg powder 3x3kg	9kg	\$275.90
CQ1000	Quantumclean laboratory detergent	1L	\$17.11
CQ1050	Quantumclean laboratory detergent	5L	\$44.75
CD1240	Decon 90 detergent	5L	\$88.65
CE2104	Extran MA 05 phosphate free	5L	\$53.40

STAINING BOXES



Used for staining, destaining and handling of electrophoresis gels and membranes. Manufactured from high quality polypropylene. Lower drain plug allows for easy liquid drainage with minimal gel



ROWE CODE	DESCRIPTION	UOM	PRICE
PS0297	Staining box 125 L x 125 W x 50mm H, for electrophoresis gels and membranes, polypropylene cover and base, LDPE plug	EA	\$8.10
PS0298	Staining box 225 L x 225 W x 50mm H, for electrophoresis gels and membranes, polypropylene cover and base, LDPE plug	EA	\$11.95

OVERSTOCKED!

K1015 Microtube high speed centrifuge.



1 ONLY ONLY AT
\$1550ea
Normally \$2225.00
SAVE \$675
IC0008

PLUS FREE ROTOR



Includes one free IR4110 Rotor valued at \$910.

Only while stocks last!

FEATURES INCLUDE:

- Accepts 0.2 - 2.2mL micro tubes
- 500 -15,000 rpm (1 rpm steps)
- RCF max of 22,000 x g
- Powerful 172 watt motor
- 0 - 99 minute timer (1 second steps)
- 235 x 235 x 350mm dimensions
- 10 inbuilt acceleration and deceleration programs
- 108 program memory

CENTURION

Limited Stock. Be Quick!

How to order online at your convenience using the Rowe Scientific Web Store and its easy to use features.



As a guest, when you arrive at the Rowe Scientific Web Store, you can search for products and populate a trolley. To obtain pricing, please request a quote on your trolley of products. Clients with a login account will see pricing directly, and will be able to place orders over the Rowe Scientific Web Store.

How to Log In

If you have a login account, enter your username and password in the login panel at the top right hand corner. This will log you into your personalised Web Store account and show your special pricing.

Your Welcome Panel

When you have successfully logged in, you will see your welcome panel to the right. Here you will find some useful features that are explained on this page.

My Backorders

My backorders is a searchable list of your current back orders with Rowe Scientific. This data is updated overnight to keep you fully informed.

Quick Order


YOUR TEN MOST PURCHASED ONLINE ITEMS

Enter quantity on any of the items in the below lists then hit the 'Add all to Trolley' button at the bottom of the page.

Part No	Description	Qty	UOM
IC3548	CONDUCTIVITY SOLUTION (0.01M KCl) (1.413e5cm @ 25C) (R934+)	50.00	1L
IC3275	CONDUCTIVITY SOLUTION (0.01M KCl) (56.676e5cm @ 25C) (R954+)	500.00	1L
FF4135	FILTER PAPER GLASS FIBRE 4.71m WHATMAN 617C (3822087)	550.00	100

Quick Order allows you place orders from a list of routinely ordered products. There are two parts to your quick order page. Your top ten recent purchases and a manual list made up of items you have added to your quick order list. These are grouped by product type e.g. Glassware for glass beakers etc.

To add a product to your quick order list, click on the add to quicklist button on the product listing. If the button does not appear on an item, it is already in your quick order list.



To order items directly by quick order list, simply enter the quantity required for each item being ordered then select the add to trolley button. You will see these items appear in your 'trolley' panel and the total value of your order.

You can continue searching and shopping for products or view your trolley and process your order.



Order History

Here you can see a list of recent online orders placed via our Web Store.



The Search Bar

The search bar functionality on the right allows you to string key words together as you would a Google search. e.g. 'beaker 250ml' the results will list those items found which best match the exact criteria to the top of the results page. When you use the product filter buttons a second inner search bar appears in the search results page allowing you to drill down from the results using key words.

Search by YOUR part number

We can upload your part numbers (where supplied by you) for contract lines allowing you to also search by your part number for easy and accurate ordering of what you need. Your part number will appear in red in brackets directly after the Rowe code.



See your agreed pricing

You will see your agreed contract and special pricing. This is indicated by the words **Your contract/special price** above the price.



Availability

You can check the availability and estimated lead times on products. Faster lead times may be available on request. The national stock is the total of our stock across the country.



Other ranges online which may interest you ...

Miele Laboratory Washers



Fever Detection Kit



Ansell Gloves



Or explore more at www.rowe.com.au

Account manager use.

USERNAME: _____ PASSWORD: _____

If you don't have a Rowe Scientific login, go to www.rowe.com.au and click the register tab to request one.



www.rowe.com.au

Online, by phone or face to face, we give you the choice.



New South Wales & ACT

Ph: (02) 9603 1205
rowensw@rowe.com.au

Queensland

Ph: (07) 3376 9411
roweqld@rowe.com.au

South Australia & NT

Ph: (08) 8186 0523
rowesa@rowe.com.au

Victoria & Tasmania

Ph: (03) 9701 7077
rowevic@rowe.com.au

Western Australia

Ph: (08) 9302 1911
rowewa@rowe.com.au



*A 100% Australian
owned company
since 1987.*

REF:547 Prices do not include GST. We reserve the right to change specifications, details and descriptions without notice. Pictures for illustrative purposes only. Discounts do not apply to service, freight and or repairs. Discount is off list price, does not apply to existing discounted pricing.
X:\MARKETING\ Lablikes.Issue8

