



A 100% Australian owned company supplying scientific laboratories throughout Australia and overseas since 1987.

LABLIKES

23

OCTOBER 2025



Ansell

MICRO-TOUCH GLOVES	SIZE	PACK SIZE	ROWE CODE
NITRILE BLUE GLOVES • Medical Examination • Powder Free • Latex Free • Non Sterile • Textured Fingers • Standard Cuff	XS	10 packs of 200 gloves	HG0767
	Small		HG0768
	Medium		HG0769
	Large		HG0770
	XL		HG0771

SPECIAL OFFER AS ANSELL ROLLS OUT THEIR INDIGENOUS SUPPORT PROGRAM
ONLY \$150
PER CARTON OF 10 PACKS OF 200 GLOVES



AUSTRALIANS in SCIENCE



Prof Richard Robson

The 2025 Nobel prize in chemistry has been awarded to professor – Richard Robson for the development of metal–organic frameworks: molecular structures that have large spaces within them, capable of capturing and storing gases and other chemicals.

Robson first discovered the metal–organic frameworks, known as MOFs for short, in 1989, with his close collaborator Bernard Hoskins.

Source: The Conversation



- ✓ One stop shop
- ✓ 100% Australian
- ✓ Online ordering
- ✓ Phone ordering
- ✓ Account Managers
- ✓ Leading brands
- ✓ Chemicals
- ✓ Consumables
- ✓ Equipment



DISCOVER LEADING BRANDS AT ROWE SCIENTIFIC PTY LTD.



Offer pricing expires 31 Jan 2026 or until stocks last and is subject to GST.

www.rowe.com.au

Page 2
Velp Advance Hotplate
Magnetic Stirrers

Page 6
Miele Glassware
Washers Bundle Deals

Velp Hot Plate Magnetic Stirrers

ADVANCE SERIES 5, 7 & 10

Our most advanced hotplate stirrer - EVER.

With all of the standard features you expect from VELP

Torque Compensation
Timer
Auto-Reverse
Intermittent Mode
Settable Rotation
ERMES Enabled
But now with: Ramps

✓ YES (speed control)
✓ YES (1 min - 99 hr 59 min)
✓ YES (5 sec - 99 hr 59 min)
✓ YES (5 sec - 99 hr 59 min)
✓ YES (CW and CCW) (control with your phone)
✓ YES
✓ UP TO 9 TEMPS ✓ UP TO 9 SPEEDS



AREX 5 ADVANCE

ROWE CODE IM1374

- 135mm diameter ceramic coated alloy plate
- Up to 20L capacity
- Heats up to 310°C



AREC 7 ADVANCE

ROWE CODE IM1375

- 180 x 180mm square ceramic coated alloy plate
- Up to 20L capacity
- Heats up to 550°C



AREX 10 ADVANCE

ROWE CODE IM1376

- 260 x 260mm square ceramic plate
- Up to 25L capacity
- Heats up to 550°C

YOU CAN ADD OPTIONS

FOR ACCURATE TEMPERATURE CONTROL ($\pm 1^{\circ}\text{C}$)

ADD PT100 TEMPERATURE PROBE ROWE CODE IP2311

WITH SUPPORT ROD ROWE CODE IR4000

AND CLAMP IF REQUIRED ROWE CODE IV0408

FOR PRECISE TEMPERATURE CONTROL ($\pm 0.5^{\circ}\text{C}$)

ADD VTF DIGITAL THERMO REGULATOR ROWE CODE IT2485

AND ROD ROWE CODE IR4000



ROWE CODES
IP2311, IR4000, IV0408

ROWE CODE
IT2485





Losing time making calibration standards?... **We have the SOLUTIONS for you!**

The purity and accuracy of standard calibration solutions are imperative to achieving not only the correct answers, but also to getting the job done efficiently, allowing you to get on with more important tasks.

Rowe Scientific are experienced in laboratory practices and understand the importance in having reliable standards and reliable supply.

SOME SOLUTIONS WE MANUFACTURE:

- ✓ AAS standards 1000ppm Ni, Fe, Cr, Mg, Co, Cu, Zn, Mn
- ✓ GOLD 2.00ppm solution in 0.1% NaOH/0.01% NaCN
ie, we can make to your matrix
- ✓ KCl conductivity standards
- ✓ Corrosion inhibitors
- ✓ All pH buffers
- ✓ Titration standards solutions: 0.500M HCl,
0.500M NaOH, 0.500M H₂SO₄, 0.1N AgNO₃
- ✓ Soil analysis buffers
- ✓ Mixed mineral acids
- ✓ Engine coolants (ethylene glycol based)

ROWE OFFER YOU:

- ✓ Primary chemicals traceable to N.I.S.T (N.B.S) standards. That's guaranteed accuracy and quality control.
- ✓ Quality Assurance Certificates with AAS standards. Available on request for other solutions.
- ✓ 100mL to 2000L quantities available.
- ✓ Custom made to meet your needs.
- ✓ Competitive pricing designed to save you money in the long term.

**Contact your local Rowe Scientific Account Manager
to find out how we can save you time and money
with Rowe chemical solution standards.**

**ACCURACY
ASSURED.
MADE BY
CHARTERED
CHEMISTS.**

Automatic Test Strip Reader QUANTOFIX RELAX

Reflectometer for evaluation of test strips

The QUANTOFIX Relax is the ideal device for the objective evaluation of MN QUANTOFIX test strips. It combines the simplicity of test strips with the precision of instrumental analysis and thus the best out of these two worlds. The QUANTOFIX Relax does not require any special strips, but evaluates the normal pH-Fix and QUANTOFIX test strips, you are currently using.

Data Memory: 200 measurements

Strips Per Hour: 50

ROWE CODE IR0370



MN MACHEREY-NAGEL

LEVO PIPETTE CONTROLLER

DLAB

0.1 - 100mL - LIMITED STOCK!

ONLY \$20 WAS \$30 SAVE \$10

- Precise Control - Drop by drop drainage
- Discharge Button - Add the last drop of liquid
- Compatible with glass or plastic serological pipettes (0.1 - 100mL)
- 3µm replaceable hydrophobic filter
- Fast Suction - Can draw 25mL of liquid in 4 seconds
- Five Colour Designs - Easy to distinguish between different applications

RED ROWE CODE DP0132

GREEN ROWE CODE DP0133

BLUE ROWE CODE DP0900

YELLOW ROWE CODE DP0134

DARK BLUE ROWE CODE DP0099




Concerned that the gasket is contaminating your analysis? You need **Chirana** 2-part disposable syringes

According to EN, ISO standards | Quality management certified according to EN ISO 13485:2016 | CE marked

- ✓ Ergonomic design, smooth flow, highly transparent barrel, tight fit
- ✓ Black extended graduation, easy to read, double stop lines, guaranteed 0mL position
- ✓ Designed tip for best connection, double stop lines, scratch-free for glove manipulation
- ✓ Pyrogen-free, latex-free, DEHP-free, PVC-free, non-toxic

VOLUME	2 PART SYRINGES	PACKING	ROWE CODE
2mL	• Luer Slip • PP Barrel • PE Plunger • Single Use (gasket free)	Bulk packed, non-sterile	PS0416
5mL			PS0417
10mL			PS0418
20mL			PS0419
2mL		Individually packed, sterile	PS0095
5mL			PS1766
10mL			PS1862
20mL			PS0359



REDIWIPE CLASSIC WHITE TOWEL

300 x 330mm (100 wipes) **ROWE CODE HT0456**

MULTI BUY OFFER

BUY 1 - 4 BOXES @ \$9.95ea

BUY 5 - 8 BOXES @ \$8.75ea

BUY 16 BOXES (2 CARTONS) @ \$8.25ea





CO₂ INCUBATOR



ROWE CODE II0161

Ambient +3°C to 55°C, CO₂ 0-20%,
humidity off or >90% RH,
IoT connectivity, touch screen
168L, includes 3 shelves

Haier

**CEO
SPECIAL
ONLY
\$7700***

Uniform and Stable Temperature:

Temperature is precisely controlled within the fluctuation range of $\pm 0.1^{\circ}\text{C}$ to ensure the normal growth of cells throughout their life cycle.

Precise CO₂ Concentration:

A new type of IR sensor controls CO₂ concentration sensitively to keep the culture medium stable in a weak alkaline environment.

180°C Dry-heat Sterilization:

Easy and effective dry heat sterilisation at up to 180°C for microorganisms including bacteria, fungi and microplasma with strong resistance in just 12 hours.

IoT APP Available for Real Time Monitoring (optional):

The operation of the incubator can be monitored at anytime through mobile APP or computer terminal, also with alarm and push function.

*While stocks lasts

UltraEco ULT FREEZER

Super energy-saving technology and precise temperature control. Designed to provide long-term storage for various biological samples and sensitive items to ensure the integrity and stability of precious materials. With its outstanding energy-saving features, it not only reduces electricity consumption but also contributes to a greener environment.

ROWE CODE	730L IF0693	830L IF0694
Temperature Uniformity ($^{\circ}\text{C}$)	± 2	± 3
ENERGY STAR Power Consumption (kWh/day)	4.16	5.18
1 - Minute Open Door Recovery Time (20°C Ambient Running at -75°C)(mins)	15	32
Warm-Up Time at 20°C Ambient (from -80°C to -50°C)(mins)	348	378
Re-Entry Time (Seconds)	60	60

49% decrease in energy consumption compared to previous generation.

**FIRST 2 UNITS ORDERED, RECEIVE A
COMPLEMENTARY SET OF 2" RACKS***
ROWE CODE IR0213

*Mention RACKS on your order



With LCD touch screen

UNDER COUNTER ULT FREEZER



ROWE CODE IF0362

-86°C, 100L, J Series

Energy Efficient, Safe and Reliable:

High efficiency HC refrigeration system, optimised for energy efficiency delivering a power consumption figure of just 5.5 kW·h/24 h

Personal ULT Storage:

810mm cabinet height makes it easy to place on or under bench, saving storage space. Stackable design.

Ergonomic Design:

Ergonomic handle design ensures easy one-hand door opening.

Low Noise:

Optimized noise reduction cabinet and system design, emits sound level of only 46.8 dB.

*While stocks lasts

**CEO
SPECIAL
ONLY
\$7250***

END OF YEAR SPECIALS

MIELE GLASSWARE WASHER BUNDLES

BUNDLE PACKAGE OPTION ONE

GLASSWARE WASHER - MODEL: PG8583 AD

ROWE CODE ID0236

PLUS EASYLOAD BASKET BUNDLE

ROWE CODE IG0055

PLUS COMMISSIONING/INSTALLATION*

OF NEW WASHER + BASIC TRAINING

ROWE CODE IC0502

Bundle offer \$14,688 SAVE \$1,100



***INSTALLATION INCLUDES 1 COMPLIMENTARY GIFT PACK OF DETERGENTS**

This offer applies to bundle offers one and two only.



BUNDLE PACKAGE OPTION TWO

GLASSWARE WASHER - MODEL: PG8583 ADP

ROWE CODE ID0104

PLUS EASYLOAD BASKET BUNDLE

ROWE CODE IG0055

PLUS COMMISSIONING/INSTALLATION*

OF NEW WASHER + BASIC TRAINING

ROWE CODE IC0502

Bundle offer \$14,845 SAVE \$1,100

BUNDLE PACKAGE OPTION THREE

GLASSWARE WASHER - MODEL: PG8504 ADP

ROWE CODE CD1543

PLUS EASYLOAD BASKET BUNDLE

ROWE CODE IG0055

PLUS COMMISSIONING/INSTALLATION

OF NEW WASHER + BASIC TRAINING

ROWE CODE IC0502

PLUS GUARDIAN CLEAN WASHING

DETERGENT TABLETS

ROWE CODE IG0157

Bundle offer \$12,000 SAVE \$800



Miele

Orders must be placed by the 17th of December 2025. While stocks last. Delivery must be completed within 60 days of order being placed. Not to be used in conjunction with any other offer. Installation and commissioning must be within 100km of the GPO on all 3 bundles offered.




DISTILLATION UNITS FOR WINE & BEVERAGE LABS

GLASSCHEM Compact distillation units for the determination of free and bound sulphur dioxide, volatile acids and alcohols are supplied pre-assembled on stainless steel mounting frames, ready to connect to your existing services and are available in different configurations to suit different laboratories' requirements. All instruments are based on accepted primary methodologies.

COMBO UNITS

OH, SO₂ and VA modules can be incorporated into "Combo" models so that multiple parameters can be measured on a single instrument to suit your laboratories' workload.



VA DISTILLATION UNIT

Available in 1, 3, 4 or 6 places



OH DISTILLATION UNIT

Available in 1, 2, 3, 4 or 6 places



SO₂ DISTILLATION UNIT

Available in 2, 4, 6 or 8 places

A full range of wine testing apparatus, alcohol stills, distillation equipment, accessories and more, can be purchased directly from your local Rowe Scientific office in Australia.

BLADEFLASK

Scalpel Blade Remover & Sharps Container


Remove, contain, and dispose of scalpel blades in one click with the BladeFLASK. The world's pioneering single-handed scalpel blade remover enables healthcare professionals to quickly and safely remove up to 100 scalpel blades.

- ✓ Puncture-proof, compliant with UN recommendations for handling dangerous waste
- ✓ Audible "click" to signify scalpel blade removal
- ✓ Built-in counter with automatic shut-off when maximum capacity is reached
- ✓ OSHA compliant with single-handed activation when used with the reusable mounting bracket
- ✓ Mount on a wall or desk with the reusable mounting bracket
- ✓ Compatible with most brands of general surgery blades and flat handles

ROWE CODE HB0500



NEW FROM BRAND GERMANY TRANSFERPETTE® PRO THE BEST PIPETTE FROM BRAND, EVER!



TRANSFERPETTE PRO ADJUSTABLE PIPETTOR

VARIABLE VOLUME	ROWE CODE
0.1 - 1µL	DT0387
0.1 - 2.5µL	DT0388
0.5 - 10µL	DT0389
2 - 20µL	DT0390
2 - 20µL	DT0391
5 - 50µL	DT0392
10 - 100µL	DT0393
20 - 200µL	DT0394
30 - 300µL	DT0395
100 - 1000µL	DT0396
250 - 2500µL	DT0397
500 - 5000µL	DT0398
1000 - 10000µL	DT0399

8 and 12 channel Transferpette PRO also now available.

TRANSFERPETTE® PRO & TRANSFERPETTE® S IN COMPARISON

	TRANSFERPETTE® PRO	TRANSFERPETTE® S
Volume range single-channel (adjustable)	12 sizes: from 0.1µL to 10mL	10 sizes: from 0.1µL to 10mL
Volume range multi-channel	6 sizes: from 0.5µL to 1250µL	5 sizes: from 0.5µL to 300µL
Size of handle	Slim	Standard
Adjustable finger rest	✓	-
One-handed operation	✓	✓
Universal tip cone	✓	✓
Volume adjustment protection	✓	✓
Counter protection	✓	-
Easy Calibration Technology	✓	✓
User Adjustment	✓	-
Completely autoclavable (121 °C)	✓	✓
UV-resistant	✓	✓
Warranty for single-channel pipettes	3 years	3 Years



INTRODUCING THE **NEW** **CHARM EZ PROTECT+ SYSTEM**

Customisable, AUTOMATED dairy antibiotic testing system.



*Charm Protect Plus has ability
to run up to 6 samples with
the addition of extra optical
modules and reader extenders.*

*Contact your local Rowe
Scientific account manager
for more details.*

PATENTED, INCUBATE AND READ IN ONE STEP!

ROSA BARCODED TESTS

Instantly ID test type, incubation time, and temperature. Identify lot# and expiration date preventing the use of expired tests.

SMART OPTICS WITH SLOPED DESIGN

Mitigate pipetting errors, identify previously run tests, and monitor proper milk flow to automatically start incubation cycle and verify proper test development.

COMPLETE DATA INTEGRATION

Integrated email to automatically notify users of various test triggers (positive results, duplicate Sample IDs, invalids, etc). Customisable, automated daily email reports.

INTELLIGENT WORKFLOWS

Users are guided through testing procedures, performance monitoring, and confirmation of positives to confidently accept or reject milk.

REAL-TIME DATA TRANSFER

Connect to LIMS, ERP, or Cloud. Connect with barcode and RFID scanners, WiFi, Bluetooth, or cellular.

ROWE CODE GC0705



Please contact your local Rowe Scientific office to arrange a complimentary in-house demonstration.

MAGNETIC HOTPLATE STIRRER

STOCK ONLY CLEARANCE

velp®

AREC

15L capacity, 550°C, 50 - 1500rpm, square ceramic top plate 180 x 180mm

ROWE CODE IM0686

STOCK ONLY CLEARANCE PRICE \$615



**6
ONLY
AT THIS
PRICE!**

AREX 6 Connect PRO

Ermes cloud platform via WiFi connection, PT100 probe, support rod and clamp, 370°C, up to 1700rpm

ROWE CODE IM0408

STOCK ONLY CLEARANCE PRICE \$1255



**8
ONLY
AT THIS
PRICE!**

AREX 6 Digital PRO

PT probe, support rod and clamp, inbuilt timer, auto reverse, 370°C, up to 1700rpm

ROWE CODE IM0557

STOCK ONLY CLEARANCE PRICE \$1095



**7
ONLY
AT THIS
PRICE!**

AREX-6 Digital PRO

370°C (0.5°C accuracy), up to 1700rpm, VTF, support rod and speed servo control

ROWE CODE IM0405

STOCK ONLY CLEARANCE PRICE \$1425



**1
ONLY
AT THIS
PRICE!**

AREC.X 7

15L capacity, 550°C, 30 - 1500rpm, square ceramic top, with probe, rod and clamp

ROWE CODE IM0696

STOCK ONLY CLEARANCE PRICE \$810



**2
ONLY
AT THIS
PRICE!**

AREC.X 7

Square ceramic top, up to 550°C, 30 - 1500rpm, socket for external probe, 15L stirring capacity

ROWE CODE IM0695

STOCK ONLY CLEARANCE PRICE \$895



**8
ONLY
AT THIS
PRICE!**



A 100% Australian owned
company since 1987.

ISSUE 23. OCTOBER 2025

ROWE SCIENTIFIC PTY LTD


LABLIKES

TEST TUBE RACK

5 x 12 position (60 place) for 16mm tubes,
no-wire, white polypropylene

ROWE CODE PR1150

**STOCK ONLY
CLEARANCE
PRICE \$34.95 EA**



LIMITED
STOCK
AT THIS
PRICE

ERLENMEYER FLASK

2000mL, wide neck, heavy duty,
heavy beaded rim, borosilicate glass

ROWE CODE GF0541

**STOCK ONLY CLEARANCE
PRICE \$19**



LIMITED
STOCK
AT THIS
PRICE

SEPARATING SQUIBB FUNNEL 250mL

With PTFE Stopcock

ROWE CODE GF4607

**STOCK ONLY CLEARANCE
PRICE \$10 EA**



LIMITED
STOCK
AT THIS
PRICE

LABORATORY BOTTLE 250mL

Amber stained GL45 borosilicate glass
(without screw cap)

ROWE CODE GB8007

**STOCK ONLY CLEARANCE
PRICE \$19 EA**



MEASURING CYLINDER 25mL

'A' class, tall form, spouted,
glass with hexagonal base

ROWE CODE GC1064

**STOCK ONLY
CLEARANCE PRICE \$10 EA**



5
ONLY
AT THIS
PRICE!



Why not order online
at your convenience
24/7 through the
Rowe Scientific Web
Store, with its easy to
use features?

Rowe Scientific Pty Ltd are a 100% Australian owned company which have been supplying laboratories across Australia since 1987 with scientific equipment, consumables and chemicals. You have the option of online ordering at your convenience 24 /7 through a secure login account we create for you. Once you have a Rowe Scientific Pty Ltd account, simply go to www.rowe.com.au log-in and then you can:

- Order online with your secure login.
- See special pricing and availability.
- Order when it's convenient for you 24 /7.
- See your online order history and back orders.
- Shop from over 117,000 product lines and leading brands.



To obtain a Rowe Scientific login, go to www.rowe.com.au and click the register tab then follow the prompts.



To view our webstore as a guest go to **www.rowe.com.au** and explore the products tab.

To see our current special offers on our webstore, go to www.rowe.com.au and explore the specials tab.

Or to find out more about us including scholarships and news events, go to www.rowe.com.au and explore the news tab.



For all the laboratory
supply needs that
you envision.



Online, by phone or face to face, we give you the choice.

INSTRUMENTS - CONSUMABLES - CHEMICALS - SERVICE & REPAIRS

A 100% Australian owned company, supplying scientific laboratories since 1987.

Discover more at: **www.rowe.com.au**

South Australia & NT
Ph: (08) 8186 0523
rowesa@rowe.com.au

Queensland
Ph: (07) 3376 9411
roweqld@rowe.com.au

Victoria & Tasmania
Ph: (03) 9701 7077
rowevic@rowe.com.au

New South Wales
Ph: (02) 9603 1205
rowensw@rowe.com.au

Western Australia
Ph: (08) 9302 1911
rowewa@rowe.com.au

