



ROWE SCIENTIFIC PTY LTD www.rowe.com.au

For accuracy and professionalism

Providing laboratory supplies to the scientific community across Australia since 1987.

We are proudly a 100% Australian owned company.



Is the chair you're sitting on up to standard? If not then...

YOU MAY WANT TO SIT DOWN FOR THIS



We're introducing a line of flexible and ergonomic chairs, **specifically made for laboratories.**



LABSIT LABORATORY CHAIR



Made for laboratories

Labsits' hygienic design makes it perfect for multiple sectors, from pharmaceutical, biotech, chemical, health care, clean room and medical technology fields.

Materials matter

In the laboratory, the perfect cleaning process starts with the material. Labsits' back rest is made from faux leather with 3D flex function. The five star base is available in polished aluminium or black plastic.

Anti bacterial and easy to clean

Made from faux leather padding which is highly effective against germs. Hygienic design with minimal joints for simple cleaning and disinfecting.

www.rowe.com.au

We have product specialists ready to help you choose your next laboratory chair



NEON LABORATORY CHAIR



Comfort for the laboratory

Neon supports the forward-leaning sitting position that is often necessary when working in the laboratory.

Materials matter

Neon is made of solid steel with an aluminium base. Also available with a choice of castors with load-sensitive brakes for hard floors or with abrasion-resistant glides.


Neon Integral foam upholstery

Extremely robust, durable, washable, resistant to mechanical damage, mild acids and alkalis.



Interchange your upholstery with ease

The Neon 1+1 System



South Australia & NT
Ph: (08) 8186 0523
rowesa@rowe.com.au

Queensland
Ph: (07) 3376 9411
roweqld@rowe.com.au

Victoria & Tasmania
Ph: (03) 9701 7077
rowevic@rowe.com.au

New South Wales
Ph: (02) 9603 1205
rowensw@rowe.com.au

Western Australia
Ph: (08) 9302 1911
rowewa@rowe.com.au