



ROWE SCIENTIFIC
PTY LTD www.rowe.com.au

For Accuracy & Professionalism

Scientific supplies, Instrument sales, Service and repair, Consulting

A proudly 100%
Australian owned
company.

ABN 63 009 437 790



CERTIFICATE OF SPECIFICATION

REAGENT:	QUANTUMCLEAN LABORATORY DETERGENT
PRODUCT CODES:	CQ1000, CQ1050, CQ1055, CQ1100, CQ1110
DESCRIPTION:	Cleaning agent concentrate for laboratory and industrial use.
APPLICATION:	For total cleaning of all laboratory and industrial glassware, plastics, ceramics and stainless steel. Not suitable for use on aluminium surfaces. Tested for decontamination of radionuclides from glassware.
APPEARANCE:	Clear to cloudy white non-viscous liquid.
FORMULATION:	Quantum Clean is a concentrated mixture of surfactants, chelating agents, stabilising agents and alkali salts in aqueous solution. The mixture is non-flammable, completely water miscible and biodegradable. Quantum clean does not contain phosphates or poly-fluoroalkyl substances.

Parameter	Specification
Density	1.02 -1.03 @ 20°C
Chelating agent content (complexometric, Zn)	>0.02 N
Free alkali (as NaOH)	1.9 - 2.1 %
Phosphate	<0.002 %

TRANSPORT INFORMATION:	Not classified as dangerous goods according to ADG7 criteria.
-------------------------------	---

NEW SOUTH WALES & ACT

Unit 1, 12 Essex Street
Minto, Sydney
New South Wales 2566
Tel: 02 9603 1205
Fax: 02 9603 1084
rowensw@rowe.com.au

QUEENSLAND

55 Westgate Street
Wacol, Brisbane
Queensland 4076
Tel: 07 3376 9411
Fax: 07 3376 9455
roweqld@rowe.com.au

SOUTH AUSTRALIA & NT

Lot 7, 41-47 O'Sullivan
Beach Rd, Lonsdale,
South Australia 5160
Tel: 08 8186 0523
Fax: 08 8186 0524
rowesa@rowe.com.au

VICTORIA & TASMANIA

Unit 2, 42 Green Street
Doveton, Melbourne
Victoria 3177
Tel: 03 9701 7077
Fax: 03 9769 2560
rowevic@rowe.com.au

WESTERN AUSTRALIA

11 Challenge Bvd
Wangara, Perth
Western Australia 6065
Tel: 08 9302 1911
Fax: 08 9302 1905
rowewa@rowe.com.au

Supplying scientific laboratories throughout Australia and abroad since 1987.